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Made in China - Fabricado na China
Fabricado en China - Προϊόντα στην Κίνα

INSTRUCTIONS FOR USE

PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USE AND KEEP THEM FOR FUTURE REFERENCE.

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Thank you for choosing Pic Solutions mobileRAPID, an automatic blood pressure monitor designed and manufactured according to the most up-to-date technologies. Before using the device, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

INTRODUCTION


The regular use of a Blood Pressure Monitor allows you to monitor blood pressure continuously, but it should not take the place of regular check-ups by your doctor. We recommend that you visit your doctor regularly for general check-ups and for more detailed information on blood pressure. Blood pressure measurements determined with this unit are equivalent to those obtained by a trained observer using cuff/ stethoscope auscultation method, within the limits prescribed by the international ESH Validation protocol.

WARRANTY

The product is guaranteed for 5 years against all conformity defects under the normal utilisation conditions provided in the user instructions. Therefore, the warranty shall not apply in case of damages caused by improper use, wear, or accidental events. Batteries, cuff and wearing parts are not included in the warranty. In case of malfunction, contact exclusively the authorised Pic Technical Customer Service.

RAPID-TECH Comfort Experience™

This unit uses the oscillometric method to detect your blood pressure. RAPID-TECH Comfort Experience™ is the latest generation technology chosen by Pic for the blood pressure monitors of the line RAPID. It allows fast and, therefore, less annoying monitoring because the measurement is performed already during the armband inflation phase. Pic mobileRAPID has been clinically validated in accordance with the ESH (European Society of Hypertension) protocol, as a guarantee of its accuracy and precision.

List of included items: Monitor with the M-L size cuff, mains adaptor, user manuals and carrying bag.

Features:

- 24,5x86,5 mm Blue LCD display with blue backlight
- Two memory areas
- 100 recorded measurement per area
- Inflation-type measurement
- Wireless connectivity (Bluetooth® version v. 4.0)
- App available for iOS and Android

ABOUT BLOOD PRESSURE

What is blood pressure? Blood pressure is the pressure exerted on the artery tube while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is diastolic (lowest) blood pressure.

Why measure your blood pressure?


Among the various health problems afflicting modern people, problems associated with high blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressures a necessity of identifying those at risk.

Blood Pressure Fluctuates

Blood pressure fluctuates all the time! You should not be over worried if you encountered two or three measurements at high levels. Blood pressure changes over the month and even throughout the day. It is also influenced by the season, temperature and activities.

Classification of high blood pressure

The out-of-office blood pressure classification published by ESH/ESC in 2013 is as follows:

Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the scheme: the scheme is meant only to discriminate between the different levels of blood pressure.

CAUTION

Only a physician can tell your normal BP range. Please contact a physician if the average of your measuring results falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector

This blood pressure monitor is equipped with an algorithm for Irregular Heartbeat (IHB) Detection. During each measurement, this algorithm records the heartbeat intervals and works out the standard deviation. If the calculated value is larger than or equal to 15, this algorithm will light up the IHB symbol on the screen when displaying the measuring result.

CAUTION

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. We wish to underline that the device does not replace a visit to the cardiologist, but it is useful to detect any irregularities in advance.

Why the blood pressure I get from the hospital is different from home?


The blood pressure is different even during 24 hour because of the weather, emotion, exercise etc, specially the "white coat" in hospital which makes the results are higher than the ones at home.

If the result is the same if measuring on the right arm?

Measurements can be taken on either arm; results can however differ sometimes so it is always advisable to use the same arm.

LCD display explanation

At starting check the completeness of the LCD segments


SYMBOL	DESCRIPTION	EXPLANATION
SYS	Systolic pressure	Maximum/systolic pressure result
DIA	Diastolic pressure	Minimum/diastolic pressure result
Pul/min	Pulse	Pulse/minute
▼	Deflating	The arm cuff is deflating
mmHg	mmHg	Blood pressure unit of measurement mmHg
⦿ + Lo	Low battery	The battery is low and needs recharging
♥	Irregular heartbeat	Irregular heartbeat
↑	Grade	Blood pressure level indicator
00:00	Date	"M" shows the month, "D" shows the day
LAST 3 AVG.	Average	The average of last 3 measurements taken within 30 minutes.
📶	Data transmitting	Data transmission succeeds (Bluetooth® functional icon)
📄	Memory	If 📄 shows, the displayed measurement values is from the memory
♥	Heartbeat	Heartbeat Detection during the measurement
👤	Motion detected while measuring	Measure could be inaccurate

The steps that the patient/user/operator can do safely are only those specified in the following instructions of use, maintenance and warnings.

Recharge the device battery

1. The battery of the mobileRAPID is a built-in rechargeable Li-polymer battery type, the battery is 3.7V 1000 mAh.
2. Please use only the AC/DC adaptor Model: KH0601000BW-PP-40 REF. 0200913500000(UK) or model KH0601000EW-PP-40 REF. 0200913600000(EU) to charge the battery, just like the following picture:

Charging the power under following circumstances:

- **⦿ + Lo** is displayed ,
- The LCD display dims,
- When powering on the monitor, the display does not light up.

CAUTION!

WARNING! Use only the PIC's authorized AC/DC Adaptor (KH0601000BW-PP-40 REF. 0200913500000(UK) or model KH0601000EW-PP-40REF. 0200913600000(EU)) to charge the power. You cannot use the blood pressure monitor during the process of charging.


WARNING! mobileRAPID has a built-in rechargeable Li-polymer battery, and cannot be replaced or disassembled by the customer and unauthorized maintenance personnel.

- The battery must be disassembled and replaced only by PIC authorized maintenance personnel. Do not try to change the battery by yourself; improper substitution or the use of wrong spare part can cause risk of explosion.
- Under the normal using, it can charge power about 300 times, if the battery cannot charge the power normally or the blood pressure monitor cannot use normally, please connect with the authorized maintenance personnel. If measured three times per day, and the battery is fully charged, it can be used for about 40 days.
- In order to extend the battery charge, please switch off the device after usage.
- Storage and use the blood pressure monitor at the cool, dry and ventilated environment.
- **Caution!** Avoid to approach to the fire and the heat source, or it will cause the battery explode. The battery shall be disposed according to the guideline shared at the end of this manual (2006/66/EC)
- During the process of charging, the blood pressure monitor displays :
 - When the charging is finished, please pull the plug.
 - When charging, shall not touch charging connector and the patient simultaneously.
- When charging with adaptor, the socket-outlet shall be installed near the equipment and shall be easily accessible.
- Lithium batteries require charging prior to use. Always refer to the instructions for correct charging procedures.
- Do not keep charging the lithium battery if fully charged.

Setting Date and Time

To ensure the stored measurement result has correct time record, please set time and unit before device is used. Before use, switch the ON/OFF button to the "I" side to turn on the monitor. Note: If the button is on the "O" side, there is no reaction when you press any button.

- (1) When the monitor is off, press and hold User 1 button for 3s to enter Time Setting Mode.
 - (2) The year will flash on the monitor, press "Query" button to change the number. Each press will increase the numeral by one in a cycling manner.
 - (3) Press "User 1" button again to confirm [YEAR]. Then the number representing the [MONTH] blinks. Press "Query" button to change the number. Each press will increase the numeral by one in a cycling manner.
 - (4) Press "User 1" button again to confirm [DAY]. Then the number representing [24 or 12h mode] blinks. Press "Query" button to change the number. Each press will change the numeral by 12 to 24h in a cycling manner.
 - (5) Press "User 1" button again to confirm [24 or 12h mode]. Then the number representing [HOURS] blinks. Press "Query" button to change the number. Each press will increase the numeral by one in a cycling manner.
 - (6) Press "User 1" button again to confirm [HOURS]. Then the number representing [MINUTE] blinks. Press "Query" button to change the number. Each press will increase the numeral by one in a cycling manner.
 - (7) Press "User 1" button again to confirm [MINUTE]. Then the LCD will display "doNE" → YEAR → DAY & MONTH → 24 or 12h → MCD → HOUR & MINUTE. Then the monitor will shut OFF.
- Note: If during the setting mode no button is pressed for more of about 60s the display turned OFF automatically.

Pairing up the Blood Pressure Monitor with Your Device

- (1) Turn on Bluetooth® connection and the app Pic Health Station. Make sure both are on when pair-up is proceeding.
- (2) When the monitor is off, press and hold the User 2 button to start pair-up. The symbol **📶** and the symbol **📄** will be shown on the LCD alternatively, indicating pair-up is proceeding. If succeed, symbol **📶** will be shown on the LCD. If fail, symbol **📄** will be shown on the LCD.
- (3) The monitor will shut off automatically after pair-up process is complete.
- (4) During the pairing, only the smartphone to be connected with the device shall run PIC APP. If you meet problem in pairing please check that no other smartphone in the family are running the APP, or reboot your smartphone.

Bluetooth® Module No. : nRF51822

Frequency Range: 2.402 - 2.480 GHz
Supply Voltage: 3.3V
Output Power Range: 4 dBm

Applying the Cuff

1. Remove all watches, jewellery, etc. prior to attaching the arm monitor. Clothing sleeves should be rolled up and the cuff should be wrapped on bare skin for correct measurements.
2. Apply cuff preferably to left arm with palm facing up as Fig. A.
3. Make sure the edge of the cuff is about 2-3 cm from the elbow.
4. In order to ensure accurate measurements, fasten the velcro strap securely around your arm so there is no extra space between the cuff and your skin. If the cuff is too loose, the measurement will not be accurate.


Fig. A

Helpful Tips


Here are a few helpful tips to help you obtain more accurate readings:

- Do not measure your blood pressure immediately after consuming a large meal.
- To obtain more accurate readings, please wait one hour after meal before measuring.
- Do not smoke or drink tea, coffee or alcohol before measuring your blood pressure.
- Do not measure your blood pressure immediately after taking a bath: wait at least 20 minutes
- You should not be physically tired or exhausted while taking a measurement.
- Do not measure your blood pressure when you want to discharge urine
- It is important that you relax during measurement.
- Do not take measurements if you are under stress or tense.
- Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
- Rest for 5 minutes before measuring.
- Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
- Take the measurement in a silent room
- Do not move and do not talk when taking a measurement
- The cuff should maintain at the same level as the right atrium of the heart (Fig. A)
- Do not cross your legs and keep your feet on the ground.
- Keep your back against the backrest of the chair
- For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician

Start the Measurement

If you want to store the measurement result on User1's memory area, press the User1 button: it will automatically turn on, finish the whole measurement and then save the result.
If you want to store the measurement result on User2's memory area, press the User2 button: it will automatically turn on, finish the whole measurement and then save the result.

1. When the monitor is off, press **▲** or **▶** button to turn on the monitor.

2. Press the User1 or User 2 button to power off, otherwise it will turn off within 1 minute. If the **📶** or **♥** icon appears, please read carefully the Helpful Tips section and repeat the measurement. Notes: To stop the measurement at any time, press the **▲** or **▶** button under the display. The cuff will quickly deflate.

CAUTION!

1. When PIC mobileRAPID successfully pair-up with your SMARTPHONE, the measurement data will be automatically transmitted to it via Bluetooth® wireless connection.

WARNING!

Data transmitted and displayed on the terminal device cannot be used for direct diagnostic purposes. Please note that the blood pressure data transmitted to any smartphone (or tablet, etc.) is not intended to replace diagnosis or to give medical advices.

Recall the measurements

1. When the Blood pressure monitor is off, press the "MEM" button to see the average of the last 3 recordings taken within 30 minutes. If the recorded measurements are less than 3 as well as the last 3 recordings are taken in more than 30 minutes, you will see the latest reading.
2. The LCD will display the latest measuring result of the user ID which completes the last measurement. The record number, measuring date and measuring time will be displayed alternatively.
3. Press "Query" button to rotate the history records.
4. When in the memory mode, press the User 1 button to recall the measurement history of User 1, or press the User 2 button to recall the measurement history.
5. When no history stored for the specific user in the monitor, press "Query" button and the LCD will display as pictured to the right.

CAUTION!

WARNING! The most recent recording (1) is shown first. Every time a new recording is made, the first (1) record is assigned to it. All the others are moved by one figure (e.g. 2 becomes 3, and so on), whereas the last record (100) is deleted from the list.

Delete the measurements

- (1) When under the query mode, press and hold "Query" button for 3 seconds to clear the memory.
- (2) When the LCD display "dEL ALL", press User button to confirm.
- (3) The LCD will display "dEL doNE" and then shut off.
- (4) Repeat points from 1 to 3 pressing User 2 to delete the user 2 measurements
- (5) If you wish to stop clearing the memory, you may press the other user button to turn off the monitor, or wait until the monitor shuts off.

Maintenance

In order to get the best performance from your monitor, please follow the below instructions.

- Put in a dry place and avoid the sunshine.
- Avoid touching water, clean it with a dry cloth in case.
- Avoid the intense shaking and collision.
- Avoid the dusty and unstable temperature environment.
- Using the wet clothing to remove the dirt.
- Avoid washing the cuff.

Please use the soft cloth to clean the whole unit. Don't use any abrasive or volatile cleaners.

Calibration and service

The accuracy of this monitor has been carefully tested and it is has been designed for a long life. Indeed the monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or five years of normal use. It is generally recommended to check the unit every two years to ensure correct functioning and accuracy. Please consult the Contact information in this user manual to get in touch with your local authorized Service Center. For battery or cuff replacement contact exclusively the authorised Pic Technical Customer Service.

WARNINGS

- This unit is not suitable for continuous monitoring during medical emergencies or operations.
- This device cannot be used with HF surgical equipment at the same time.
- This device is not intended for patient transports outside a healthcare facility.
- This device is intended for adult use only. Keep out of reach of children.
- People with reduced cognitive abilities shall use the device under supervision.
- This device is intended for non-invasive measuring and monitoring of arterial blood pressure.
- It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- Do not confuse self-monitoring with self-diagnosis. By using this device you can monitor your blood pressure, but it does not replace your doctor. Begin or end medical treatment based solely on medical advice. If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your Physician.
- If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a signal of will be displayed. Under this condition, the Electronic Sphygmomanometer can keep function, but the results may not be accurate, it's suggested that you consult with your physician for accurate assessment
- When the device is used for measurements on patients with common arrhythmias, such as atrial or ventricular premature beats or atrial fibrillation, a deviation may occur. Please consult your doctor about the result. If the cuff pressure exceeds 300 mmHg it deflate automatically. If the arm cuff does not deflate when pressure exceeds 300 mmHg, take it off your arm and press the **▲** or **▶** button to stop inflation.
- The device is a non-AP/APG device and not suitable for use in the presence of anaesthetic mixtures that are inflammable with air, oxygen or nitrous oxide.
- The operator shall not touch DC output jack of AC/DC adapter and the patient simultaneously.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- This device is contraindicated for any female who may be suspected of, or is pregnant. Besides provided inaccurate readings, the effects of this device on the fetus are unknown.
- Manufacturer will make available on request circuit diagrams, component parts list etc. Only to qualified person.
- Too frequent measurements can cause injury to the PATIENT due to blood flow interference. Please check (by observation of the limb concerned) that operation of mobileRAPID does not result in prolonged impairment of the circulation of the blood of the PATIENT.
- This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become anaesthetic, swollen and even purple due to a lack of blood.
- Please use the device under the environment which was provided in the instruction manual. If this is not complied with, correct operation and life of the device may be compromised.
- During using, the patient will contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of ISO 10993-5:2009 and ISO 10993-10:2010. It will not cause any potential sensitization or irritation reaction.
- The application of the CUFF over a wound, as this can cause further injury
- The application of the CUFF and its pressurization on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present because of temporary interference to blood flow and could result in injury to the PATIENT
- The application of the CUFF and its pressurization on the arm on the side of a mastectomy can cause injury or give not accurate measurement

- The pressurization of the CUFF can temporarily cause loss of function of simultaneously used monitoring ME EQUIPMENT on the same limb.
- Please use accessories and detachable parts specified/ authorised by the manufacturer. Otherwise, it may cause damage to the unit or danger to the user/patients.
- We recommend checking its performance every two years or after it has been repaired. The device must be checked again in case of damage due to impact (such as falls) or exposure to liquids and/or extreme temperatures (hot/cold) or extreme humidity variations.
- For public place use, please contact your local distributor or the Italian toll-free number 800 900 080 to get support how to proceed.
- Please dispose of accessories, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- In case of problems concerning settings, maintenance or use of the device, contact your local distributor or the Italian toll-free number **800 900 080**. Warning! Do not service and/or repair the appliance while it is being used! Any other maintenance/ service activity/operation must be done solely by a PIC Service Centre.

SAFETY INFORMATION

The symbols below can be found in the Instruction Manual as well as on labels and other components.

	"PLEASE REFER TO THE INSTRUCTIONS FOR USE" symbol		"TYPE BF APPLIED PART" symbol
	Symbol for "COMPLIES WITH MDD 93/42/EEC as amended by directive 2007/47/EC"		Storage and transport humidity symbol
	"MANUFACTURER" symbol		Off symbol
	Symbol for "SERIAL NUMBER"		On symbol
	Latex free		Caution
	Catalogue Number		The Bluetooth® combination mark, Bluetooth® version v.4.0
	Symbol for "COMPLIES WITH RED 2014/53/EU Directive"		Interference may occur in the vicinity of equipment marked with this symbol
	Simbolo symbol for "ALTERNATE CURRENT"		Symbol for "DIRECT CURRENT"
IP22	Degree of protection against penetration of solids and liquids (dripping water when tilted 15°).		

Appliance date of production: The second pair of digits in the lot number identifies the year, the third pair identifies the month. (e.g. LOT 1016100001: 16=2016; 10= October).

This section includes a list of error messages and frequently asked questions for problems you may encounter with your sphygmomanometer. If the products not operating as you think it should, check here before arranging for servicing.

PROBLEM	SYMPTOM	CHECK THIS	REMEDY
No power	Display will not light up.	Low battery. ON (I) - OFF (O) button switched OFF (O).	Recharge the batteries. switch ON (I) the ON (I) - OFF (O) button.
Low batteries	The display is dim or shows ⦿ + Lo	Low battery.	Recharge the battery.
E12 shows		Communication error	Check if the App is on or not, try data transmission again.
E01 shows		The cuff is not secure	Refasten the cuff and relax for a moment and then measure again.
E02 shows	The monitor detected motion, talking or the pulse is too poor while measuring.		Relax for a moment and then measure again.
E03 shows		The measurement process does not detect the pulse signal.	Loosen the clothing on the arm and then measure again.
E04 shows		The treatment of the measurement failed.	Relax for a moment and then measure again.
EExx shows		A calibration error occurred.	Retake the measurement. If the problem persists, contact the retailer or our customer service department for further assistance. Refer to the warranty for contact information and return instructions.
OUT		Out of range measurement.	Relax for a moment. Refasten the cuff and then measure again. If the problem persists, contact your Physician.

SPECIFICATIONS

Power supply:	3.7V 1000mAh Built-in rechargeable Li-polymer battery (~120 measurements with full battery, battery rechargeable 300 times)
Display:	Digital Backlight, LCD 25x86mm
Measurement:	Oscillographic testing mode
Measurement range:	Rated cuff pressure: 0mmHg~300mmHg
Measurement pressure:	SYS: 60-230mmHg, DIA: 40-130mmHg
Pulse value:	40-199 beat/minute
Accuracy:	Pressure: ± 3mmHg Pulse value: ±5%
Normal working condition:	Temperatures: from 5°C to 40°C Relative humidity≤85%
Storage & transportation condition:	Atmospheric pressure: 86kPa to 106kPa Temperatures: -20°C +60°C Relative humidity 10%-93%
Atmospheric Pressure:	50-106 kPa
Arm circumference:	About 22cm~42cm
Net Weight:	Approx. 295 g (including battery and cuff)
External dimensions:	Approx. 74.3x133.0x25.4mm (excluding cuff)
Mode of operation:	Continuous operation
Useful life:	10,000 Measurements
Accessories:	
AC/DC Adaptor:	100-240V~ 50/60Hz / 6V 1A AC/DC Adaptor Model KH0601000BW-PP-40 REF. 0200913500000(UK) or model KH0601000EW-PP-40 REF. 0200913600000(EU) please refer to customer service.
Degree of protection:	Type BF applied part (cuff)
Protection against ingress of water:	IP22 degree of protection against penetration of solids and liquids (dripping water when tilted 15°)
Software version:	A01
Device classification:	Battery Powered Mode: Internally Powered ME Equipment
AC/DC Adaptor:	Class II Equipment

Contact Information

For further information on our products, please visit the website www.picsolution.com. Here you can find customer service, solve ordinary problems.

EMC Guidance

- 1) mobileRAPID needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided in the ACCOMPANYING DOCUMENTS
- 2) Portable and mobile RF communications equipment can affect mobileRAPID
- 3) Wireless communications equipment such as wireless home network devices, mobile phones, cordless telephones and their base stations, walkie-talkies can affect this equipment and should be kept at least a distance d=3.3m away from the equipment. The distance d is calculated by the manufacturer from the 800 MHz to 2.5 GHz column

Guidance and manufacturer's declaration - electromagnetic emissions		
The mobileRAPID is intended for use in the electromagnetic environment specified below. The customer or the user of the mobileRAPID should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The mobileRAPID uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. mobileRAPID is suitable for use in every room or domestic environment, and in rooms directly connected to a low voltage public network powering buildings used for domestic purposes.
RF emissions CISPR 11	Class [B]	
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic emissions			
The mobileRAPID is intended for use in the electromagnetic environment specified below. The customer or the user of the mobileRAPID should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV ai	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines	Mains power quality should be that of a typical commercial or hospital environment. The electrical fast transient burst (EFT) is generated by the switching of inductive loads. Separation between the equipment and other loads shall be considered before installation. Mains filter is required, if necessary.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s)	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5s	<5% UT (>95% dip in UT) for 0.5 cycle 40% UT (60% dip in UT) for 5 cycle 70% UT (30% dip in UT) for 25 cycle <5% UT (>95% dip in UT) for 5s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the mobileRAPID requires continued operation during power mains interruptions, it is recommended the mobileRAPID be powered from an uninterruptible power supply or a battery.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE UT is the a.c. mains voltage prior to application of the test level.			

Guidance and manufacturer's declaration - electromagnetic emissions			
The mobileRAPID is intended for use in the electromagnetic environment specified below. The customer or the user of the mobileRAPID should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3Vrms 150KHz to 80MHz	3 V	Portable and mobile RF communications equipment should be used no closer to any part of the mobileRAPID, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = \left[\frac{3,5}{V1} \right] \sqrt{P}$ $d = \left[\frac{3,5}{E1} \right] \sqrt{P} \quad 80 \text{ MHz to } 800 \text{ MHz}$ $d = \left[\frac{7}{E1} \right] \sqrt{P} \quad 800 \text{ MHz to } 2,5 \text{ GHz}$ where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ⁴ should be less than the compliance level in each frequency range. ⁵ Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80MHz to	3 V/m	
NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
⁴ Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the mobileRAPID is used exceeds the applicable RF compliance level above, the mobileRAPID should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the mobileRAPID. ⁵ Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Recommended separation distances between portable and mobile RF communications equipment and the mobileRAPID			
The mobileRAPID is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the mobileRAPID can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the mobileRAPID as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz d=1.2	80MHz to 800MHz d=1.2	800MHz to 2.5GHz d=2.3
0,01	0.12	0.12	0.23
0,1	0.37	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	13
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

THE PRODUCT COMPLIES WITH:

EN 1060-1: 1995/A2: 2009 Non-invasive sphygmomanometers – Part 1: General requirements
EN 1060-3: 1997/A2: 2009 Non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems
EN 1060-4: 2004 Non-invasive sphygmomanometers – Part 4: Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

THIS PRODUCT COMPLIES WITH EU DIRECTIVE 2006/66/EC

The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for bringing batteries, at the end of their life, to the appropriate collection facilities in order to facilitate treatment and recycling. The correct disposal of your old toy will help prevent potential negative consequences for the environment and human health and favour the recycling of the materials of which the product is made. Illegal disposal of the product by the user will damage the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

DECLARATION OF EC COMPLIANCE:

PIKDARE S.r.l. hereby declares that this appliance, mobileRAPID complies with all the essential requirements and other provisions set forth by the EU Directives: 2014/53/EU (RED), 2011/65/EU (RoHS) and 93/42/EEC amended by 2007/47/EC (MDD).
A full copy of the EC declaration can be requested from PIKDARE S.r.l. Via S. Catelli 10 - 22070 Casnate con Bernate (CO) Italy - www.picsolution.com

NOTE:

- The Bluetooth® trademark and logos are registered trademarks of SIG, Inc. and any use thereof by PIKDARE is under license.
- Google Play and the Google Play logo are Google Inc. trademarks.
- Apple and the Apple logo are trademarks of Apple Inc., registered in the United States and in other countries. App Store is a service trademark of Apple Inc. registered in the United States and in other countries.

Pouch cleaning symbols

- Wash by hand in cold water,
- Do not bleach,
- Do not tumble dry,
- Line drying,
- Do not iron,
- Do not dry clean.

PIKDARE S.r.l.,
Via Saldarini Catelli 10,
22070 - Casnate con Bernate (CO) - Italy

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- Google Play e o logotipo Google Play são marcas registadas propriedade da Google Inc.
- App Store é uma marca registada de serviços da Apple Inc. registada nos Estados Unidos e em outros países.