
Infrared ear thermometer

## INSTRUCTIONS FOR USE

**PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USE AND KEEP THEM FOR FUTURE REFERENCE.**

Dear Customer, thank you for choosing our product!  
All PIC Solution thermometers are produced using the most advanced technology and undergo stringent quality control checks, in accordance with EU legislation, to guarantee the utmost quality and precision. Do not use this device for Live & Dead decision or Safety related applications. In case of doubts or increased body temperature, consult a physician.

**Intended Use:** The infrared ear thermometer measures the infrared heat generated by the eardrum and surrounding tissue to reflect patient's body temperature  
**Intended operator:** At least 11 years old (5 years intensive reading experience), no maximum.

### HOW TO TAKE A TEMPERATURE Switching on the appliance and taking readings

1. Press "ON/MEM" button (see figure 1), the thermometer is ready for use after the  icon stop flashing and two short beep sound.


2. Taking Ear Temperature:

2.1 Gently move the ear back slightly to straighten the ear canal and very carefully insert the probe into the canal. Aim towards the timpani membrane to get the most accurate reading (fig. 2).

- For children under two years of age, gently pull the ear backwards.
- For children over two years of age and adults, gently pull the ear upwards and backwards.

2.2 Press and hold the "MEASURE" button until you hear a beep sound. Remove the probe from ear and read the temperature measurement on display.

2.3 For Next Measurement: After two short beeps are heard, the  icon stops flashing to be ready for next measurement.

### PRECAUTIONS

- The thermometer tends to compensate ear temperature readings according to the room temperatures. To get the most accurate reading the thermometer must stay in the room where it will be used, at stable ambient temperature, for at least 15 minutes before operating. The person whose temperature will be taken should also be in the same room for at least 5 minutes before taking the reading. Whilst taking the temperature, make sure that the patient and thermometer are not exposed to direct heat, sunlight or air conditioning.
- Before taking the reading, check that the ear canal is clean and accessible; do not clean the ear just before using the probe as the use of hot or cold water will affect the reading.
- Do not take a person's temperature for at least thirty minutes after physical effort, a bath or after a meal.
- Always make sure the probe is clean and without damage.

### REMARKS:

1. It is recommended that you measure 3 times with the same ear. If the 3 measurements are different, select the highest temperature.
2. Before the measurement, please stay in a stable environment for 5mins and avoid the exercise, bath for 30mins.
3. To avoid the risk of cross contamination, please clean the probe according to "Care and cleaning" section after each use.
4. The different temperature varies in healthy persons: between different parts of the body can be between 0.2~1°C
5. The "Clinical Bias" is -0.5~ -0.7°C.
6. The "Limits of Agreement" is 1.02.
7. The "Repeatability" is 0.19°C.

**WARNING:** Choking from swallowing small parts and batteries by children or pets is possible, please keep small parts and batteries at places where children and pets can't reach.

### PAIRING UP THE THERMOMETER WITH YOUR SMARTPHONE

1. Download the free app PIC HEALTH STATION from the Google Play and AppStore.
  2. Turn on PIC HEALTH STATION app on your smartphone. Bluetooth® version v. 4.0 will be automatically turned on.
  3. Make sure the thermometer is powered on.
  4. Search the thermometer with PIC HEALTH STATION and pair the thermometer with your smartphone.
  5. During the pairing, only the smartphone to be connected with the thermometer shall run PIC HEALTH STATION. If you meet problem in pairing please check that no other smartphone are running PIC HEALTH STATION, or reboot your smartphone.
- Please take note that the thermometer sends temperature data only when it is connected with your smartphone.  
Every temperature reading will be transfer to your mobile device automatically.

### AFTER MEASUREMENT

1. Power off: Device will automatically shut off if left idle for more than 1 minute to extend battery life.
2. Clean the probe after each use to ensure an accurate reading and avoid cross contamination. (see section of Cleaning and Storage for details.)
3. Remember to put on the Probe Cap when not in use.

### FUNCTIONS

#### Fever Indication

If the thermometer detects a body temperature  $\geq 37.5^{\circ}\text{C}$  (99.5°F), there will be "Red" light signal. At the same time, a long beep sound followed by three short beep sounds to warn the user for potential fever; if temperature  $< 37.5^{\circ}\text{C}$ , there will be "Green" light signal.

**WARNING:** Red light - suggest continually monitoring the body temperature. If in doubt about the person's condition, seek medical assistance.

#### MEMORY

There are total 25 sets memories for body temperature.  
✳ When power on, press the "ON/MEM" button to see the temperature stored with  icon. Please take note that data stored in the memory of the thermometer are not sent to the app.

#### °C / °F SWITCH

In "Power off" mode. Press and hold the "MEASURE" button, then immediately press and release the "ON/MEM" button. After 3 seconds, icon  will be switched to icon . (You may use the same way to switch back to °C)

#### MUTE MODE

The device setting with buzzer is on, you can set buzzer on/off under Mute mode.

When power on, press and hold the "ON/MEM" button for 3 seconds. The icon  will flash on the LCD screen and then release the "ON/MEM" button to set MUTE. Thus you will not hear beep sounds. You can also use the same process to turn off the Mute function.  
**NOTE:** If keep pressing "ON/MEM" button for 5 seconds after icon flashing, the device will be power off WITHOUT setting Mute.

#### Bluetooth®

The temperature data is automatically synced to the device via Bluetooth® v. 4.0

#### SPECIFICATIONS

- Comply with ASTM E1965-98, EN ISO 80601-2-56, IEC/EN60601-1-2(EMC), IEC/EN60601-1(Safety) standards, ISO10993, RoHS
- **Temperature measurement range:** 34~42.2°C (93.2~108°F)
- **Operating range:** 10~40°C (50~104°F), 15%~85% RH
- **Storage/Transportation environment:** Stored at room temperature between -20~+50°C, RH  $\leq$  85%
- **Atmospheric pressure:** 800~1013 hPa
- **Accuracy:**  $\pm 0.2^{\circ}\text{C}$  (0.4°F) within 35~42°C (95~107.6°F),  $\pm 0.3^{\circ}\text{C}$  (0.5°F) for other range.
- **Battery:** One lithium cell battery (CR2032 \*1 pcs)
- **Enclosure Rating:** IP54
- **IP54 degree:** Protected against splashing water and dust-protected.
- **Dimensions:** 143.16 x 39.91 x 55.86 mm
- **Weight:** 75.7 grams including battery
- **Battery life:** Around 2,000 continuous readings.
- **Expected Service Life:** 4 years
- **BT module:** nRF51822
- **Frequency:** 2,402 - 2,480 GHz
- **Output power range:**  $\leq 4$  dBm
- Bluetooth® & Fever indication & memory and °C/°F switch function.
- This thermometer converts the ear temperature to display its "oral equivalent." (according to the result of the clinical evaluation to get the offset value).

#### WARNING:

- The device should not submerge into any liquids and expose it to direct moisture.
- There is no gender and age limitation for using infrared thermometer.
- This is not an AP or APG product and is therefore not suitable for use in the presence of mixes of flammable or explosive anaesthetic gases, air, oxygen or nitrous oxide..

#### CARE AND CLEANING

**WARNING:** The probe is the most delicate part of the thermometer. Use with care when cleaning the lens to avoid damage.

1. After the measurement, please use the cotton swab with the Alcohol (70% concentration) to clean the lens (on the inside of the probe).
2. Allow the probe to fully dry for at least 1 minute.
3. Storage temperature range: It should be stored at room temperature between -20~+50°C, RH  $\leq$  85%
4. Keep the unit dry and away from any liquids and direct sunlight.
5. The Probe should not be submerged into liquids.

**\*Note:** Please check the device if damaged once it falls. If you can't make sure of it, please send the complete device to the nearest retailer for recalibration.

**WARNING:** Holding the thermometer too long may cause a higher ambient temperature reading of the probe. This could make the body temperature measurement lower than usual.

#### CHANGING THE BATTERY

- Battery replacement 
- This device is supplied with one lithium cell (CR2032x1).

1. Please pick the rubber gasket on the battery cover by using small, pointed screwdriver with the "X" shaped.
2. Loosen the screw on the battery cover.
3. Use thumb to push battery cover out.
4. Hold the device and use the screwdriver to flip out the battery.
5. Replace new battery under the metal hook on the left side and press the right side of the battery down until you hear a "click" sound.
6. Replace the battery cover and fasten the screw, and stuff the screw hole with the rubber gasket.

**WARNING:** The positive (+) side up and the negative (-) side pointed down.

#### TROUBLE SHOOTING

Error Message	Problem	Solution
	Error 5~9, the system is not functioning properly.	Unload the battery, wait for 1 minute and repower it. If the message reappears, contact the retailer for service.
	Measurement before device stabilization.	Wait until all the icons stops flashing.
	The ambient temperature is not within the range between 10oC and 40oC (50F ~ 104oF).	Allow the thermometer to rest in a room for at least 15 minutes at room temperature: 10oC and 40oC (50oF ~ 104oF).
	Temperature taken is higher than +42.2oC (108°F)	Spend 5 minutes at room temperature and avoid physical exercise and bathing for 30 minutes; then make sure that the probe lens is clean before taking the temperature again.. If a malfunction still exists, please contact the nearest retailer.
	Temperature taken is lower than +34oC (93.2°F)	
	Device cannot be powered on to the ready stage.	Change with a new battery.

#### SYMBOLS KEY

	The CE mark and Notified Body Registration Numbers, the requirement of Annex II from Medical Device Directive 93/42/EEC are met.		Indicates this device is subject to the Waste Electrical and Electronic Equipment Directive 2012/19/EU in the European Union. To protect the environment, dispose of useless device at appropriate collection sites according to national or local regulations.
	Authorized representative in the European community		Caution
	Manufacturer		BF type applied part
	Battery Recycling		Paper Recycling
	Classification of enclosure protection against ingress of water and particulate matter.		The Bluetooth® combination mark, Bluetooth® version v. 4.0
	Follow instructions for use		Product catalogue reference code or catalog number
	Instructions for use		Stand-by
	Batch number or Lot Number		Compliance to RED directive 2014/53/EU

Guidance and manufacturer's declaration – electromagnetic emissions		
The thermometer intended for use in the electromagnetic environment specified below. The customer or the user of the thermometer should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The thermometer uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.  The thermometer is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Not applicable	

Guidance and manufacturer's declaration – electromagnetic immunity			
The thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the thermometer should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	Not applicable	Portable and mobile RF communications equipment should be used no closer to any part of the thermometer, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.  <b>Recommended separation distance:</b>  $d = 1,2 \sqrt{P}$  $d = 1,2 \sqrt{P}$ 80 MHz to 800 MHz  $d = 2,3 \sqrt{P}$ 800 MHz to 2,5 GHz  where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, a should be less than the compliance level in each frequency range b. Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2,5 GHz		
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the NOME PRODOTTO is used exceeds the applicable RF compliance level above, the thermometer should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the thermometer.
- b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Guidance and manufacturer's declaration – electromagnetic immunity			
The thermometer is intended for use in the electromagnetic environment specified below. The customer or the user of the thermometer should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	6 kV contact 8 kV air	6 kV contact 8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/burst IEC 61000-4-4	2 kV for power supply lines 1 kV for input/output lines	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	1 kV line(s) to line(s) 2 kV line(s) to earth	Not applicable	Mains power quality should be that of a typical commercial or hospital environment.
Interruptions and voltage variations on power supply input lines. IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0,5 cycle  40 % UT (60 % dip in UT) for 5 cycles  70 % UT (30 % dip in UT) for 25 cycles  <5 % UT (>95 % dip in UT) for 5 sec	Not applicable	Mains power quality should be that of a typical commercial or hospital environment. If the user of the thermometer requires continued operation during power mains interruptions, it is recommended that the thermometer be powered from an uninterruptible power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
NOTE: UT is the a.c. mains voltage prior to application of the test level.			

Recommended separation distances between portable and mobile RF communications equipment and the ME EQUIPMENT or ME SYSTEM			
The thermometer is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the thermometer can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the thermometer as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter M		
	150 KHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2,5 GHz $d = 2,3 \sqrt{P}$
0.01	N/A	0.12	0.23
0.1	N/A	0.38	0.73
1	N/A	1.2	2.3
10	N/A	3.8	7.3
100	N/A	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer. NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

#### DECLARATION OF EC COMPLIANCE:

Hereby, RADIANT INNOVATION INC declares that the radio equipment type REF. 02009196100000, is in compliance with Directive 2014/53/EU, Directive 93/42/EEC and directive 2011/65/EU (RoHS). The full text of the EU declaration of conformity is available at the following internet address: <http://www.radiantek.com.tw/download/certifications.asp>

#### This product complies with the Directive 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

#### COMPLIANCE WITH EU DIRECTIVE 2006/66 CE

The crossed bin symbol on the batteries or product pack indicates that, at the end of their life, they must not be disposed of as urban refuse. They must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for taking the batteries to special waste disposal sites at the end of their life to allow for their processing and recycling. The correct disposal of your old toy will help prevent potential negative consequences for the environment and human health and favour the recycling of the materials of which the product is made. Improper disposal of the product by the user may harm the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

#### NOTE:

- The Bluetooth® trademark and logos are registered trademarks of SIG. Inc. and any use thereof by PIKDARE is under license.
- Google Play and the Google Play logo are Google Inc. trademarks.
- Apple and the Apple logo are trademarks of Apple Inc., registered in the United States and in other countries. App Store is a service trademark of Apple Inc. registered in the United States and in other countries.

#### WARNING:

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions for use. The warranty shall not therefore apply in the case of damages caused by improper use, wear or accidental events.

For the duration of warranty on conformity defects please refer to the specific provisions of national laws applicable in the country of purchase, where provided.

In case of problems concerning settings, maintenance or use of the device, contact your local distributor or the Italian toll-free number 800 900 080.

**WARNING:** Do not service and/or repair the appliance while it is being used! Any other maintenance/service activity/operation must be done solely by a PIC Service Centre.

Manufacture Date: as the serial number (please open the battery cover, it is shown on the inside of the device).

LOT: 10YYMMXXXX, appliance date of production: the second pair of digits in the lot number identifies the year, the third pair identifies the month (LOT 1016100001; 16=2016; 10=October).

Note: The thermometer is calibrated at the time of manufacture. If at any time you question the accuracy of temperature measurements, please contact the your local distributor or the Italian toll-free number 800 900 080.

**WARNING:** No modification of this equipment is allowed

 Please read these instructions carefully before use and keep them for future reference

 Type BF applied part

#### WARRANTY

The product is guaranteed against any conformity defect in normal conditions of use as provided for by the instructions for use. The warranty shall not therefore apply in the case of damages caused by improper use, wear or accidental events.

For the duration of warranty on conformity defects please refer to the specific provisions of national laws applicable in the country of purchase, where provided.

In case of problems concerning settings, maintenance or use of the device, contact your local distributor or the Italian toll-free number 800 900 080. Warning! Do not service and/or repair the appliance.

Medical Technology Promed Consulting GmbH  
Add: Altenhofstrasse 80, D-66386 St. Ingbert, Germany

Radiant Innovation Inc. <http://www.radiantek.com.tw>  
Add: 1F, No.3, Industrial East 9th Road, Science-Based Industrial Park,  
HsinChu, Taiwan 300.

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Medical Technology Promedct Consulting GmbH
PIKDARE S.r.l.,