## ISTRUZIONI D'USO

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**INSTRUCTIONS FOR USE**

**Before use, carefully read the information in this manual, and keep for future reference. If there are doubts about the interpretation of any contents of this instructions manual, contact the dealer, distributor or manufacturer.**

Dear customer, Thank you for choosing Pic Solutions MobileRAPID Wrist, an automatic blood pressure monitor designed and manufactured according to the most up-to-date technologies. Before using the device, we highly recommend that you read the short list of warnings included in this manual, to make sure you have correctly understood how to use it.

**INTRODUCTION**

The regular use of a Blood Pressure Monitor allows you to monitor blood pressure continuously, but it should not take the place of regular check-ups by your doctor. We recommend that you visit your doctor regularly for general check-ups and for more detailed information on blood pressure. Blood pressure measurements determined with this unit are equivalent to those obtained by a trained observer using cuff/ stethoscope auscultation method, within the limits prescribed by the international ESH Validation protocol.

**WARRANTY**

The product is guaranteed for 5 years against all conformity defects under the normal utilisation conditions provided by the user instructions. Therefore, the warranty shall not apply in case of damages caused by improper use, wear or accidental events.

Batteries, cuff and wearing parts are not included in the warranty. In case of malfunction, contact exclusively the authorised Pic Technical Customer Service.

**RAPID-TECH Comfort Experience™**

This unit uses the oscillometric method to detect your blood pressure. RAPID-TECH Comfort Experience™ is the latest generation technology chosen by PIC for the blood pressure monitors of the line RAPID. It allows fast and, therefore, less annoying monitoring because the measurement is performed already during the cuff inflation phase. Pic MobileRAPID Wrist has been clinically validated.

**NAME/FUNCTION OF EACH PART**

**ABOUT BLOOD PRESSURE**

What is blood pressure? Blood pressure is the pressure exerted on the artery tube while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

**Why measure your blood pressure?**

Among the various health problems afflicting modern people, problems associated with high blood pressure are by far the most common. High blood pressure has dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressure a necessity to identify those at risk.


**Blood Pressure Fluctuation**

Blood pressure fluctuates all the time! You should not be over worried if you encountered two or three measurements at high levels. Blood pressure changes over the month and even throughout the day. It is also influenced by the season, temperature and activities.

**Classification of blood pressure**

The out-of-office blood pressure classification published by ESH/ESC in 2013 is as follows.

**Note:** It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the scheme: the scheme is meant only to discriminate between the different levels of blood pressure.


**WARNING**  
Only a physician can tell your normal BP range. Please contact a physician if your measuring result falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

**Why the blood pressure I get from the hospital is different from home?**  
The blood pressure is different even during 24 hours because of the weather, emotion, exercise, etc. especially the "white coat" in hospital which makes the results are higher than the ones at home.

**If the result is the same if measuring on the right arm?**  
Measurements can be taken on either arm; results can however differ sometimes so it is always advisable to use the same arm.

**LCD display explanation**

**WARNING**  
In order to protect your device and yourself and to get the best performances, use only the batteries according to the information provided in this manual (2 batteries 1.5V AAA; 3 Vdc).

**Installing and Replacing the Batteries.**

- Slide off the battery cover.
- Install the batteries by matching the correct polarity, as shown.
- Replace the cover.

**Low battery warning:**

It is necessary to replace the batteries when the Low Battery symbol appears on the display or when the display does not turn on after the POWER button is pressed.

- WARNING**
- If the device is not used for a while, please remove the batteries.
  - Old batteries are harmful to the environment, so please disposal accordingly to the indication reported in this user manual (EU Directive 2006/66/EC).
  - Do not dispose of batteries in fire. Batteries may explode or leak.
  - There are no user serviceable parts inside. Batteries or damage from old batteries are not covered by warranty.
  - Only use brand batteries. Always replace with new batteries together. Use batteries of the same brand and same type. Other types of batteries can cause explosions or fires.
  - The battery can explode if handled carelessly. Never recharge, remove or position the batteries near sources of heat. Do not incinerate.
  - Rechargeable batteries are not recommended.

**Setting the Backlight Function/Date and Time**

- It is necessary to set the date and time for the unit every time batteries are installed or replaced.
- When the unit is off press and release the SET button, "bl.01" appears on the display, with "01" blinking.
  - Press the or button to turn on or off the backlight.
    - BL 00: backlight OFF
    - BL 01: backlight ON
  - With backlight ON the monitor is always lighted during operating. The backlight turns off temporarily after 20 sec. idling to conserve power energy.
- After the backlight setting, the "Year" will begin to flash on the display.
- Press the button to increase or to decrease till reaching the desired "Year".
  - Press "START/STOP" or "Set" button to confirm the "Year".
  - Next, the "Month" will blink. Repeat steps 1 and 2 to set the "Month" and "Day", then "Hour", and "Minute".
  - After setting the minute, the unit briefly shows the word OFF before shutting down.

**Applying the Cuff**

- Remove all watches, jewelry, etc. before attaching the wrist monitor. Clothing sleeves should be rolled up and the cuff should be wrapped on bare skin for correct measurements.
  - Apply the cuff preferably to left wrist with palm facing up.
  - Make sure the edge of the cuff is about 0.5-1 cm from the palm as shown in fig. A.
  - In order to ensure accurate measurements, fasten the velcro strap securely around your wrist so there is no extra space between the cuff and the wrist.
- If the cuff is not wrapped tight enough, the measurement values may be false.

Fig. A

**Correct Measuring Posture**

- Place your elbow on the table so that the cuff is at the same level as your heart as Fig. B.
- Turn your palm upwards.
- Sit upright in a chair, and take 5-6 deep breaths.

**Helpful Tips**

- Here are a few helpful tips to help you obtain more accurate readings:
- Do not measure your blood pressure immediately after consuming a large meal.
  - To obtain more accurate readings, please wait one hour after meal before measuring.
  - Do not smoke or drink tea, coffee or alcohol before measuring your blood pressure.
  - Do not measure your blood pressure immediately after taking a bath: wait at least 20 minutes.
  - You should not be physically tired or exhausted while taking a measurement.
  - Do not measure your blood pressure when you want to discharge urine.
  - It is important that you relax during measurement.
  - Do not take measurements if you are under stress or tense.
  - Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
  - Rest for 5 minutes before measuring.
  - Wait at least 3 minutes between measurements. This allows your blood circulation to recover.
  - Take the measurement in a silent room.
  - Do not move and do not talk when taking a measurement.
  - The cuff should maintain at the same level as the right atrium of the heart (Fig. B).
  - Do not cross your legs and keep your feet on the ground.
  - Keep your back against the backrest of the chair.
  - For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

**Start the measurement**

Please read the preceding portions of this manual before taking your first reading.

- After applying the cuff and your body is in a comfortable position, press the START/STOP button to turn on the device.
- All display segments are shown for self-test. Please contact the service center if segment is missing (see for reference LCD display).
- After the self-test, the valving for the last reading appears on the display, then the monitor starts to seek zero pressure.
- The cuff automatically begins to inflate, with the display showing the increasing pressure in the cuff. As the pressure increases, an arrow pointing up appears on the display.
- When the monitor detects the heartbeat, the heartbeat symbol blinks on the LCD display.
- Finally the blood pressure and pulse rate is calculated and displayed on the LCD screen together with ESH indicator, date and time.
- Press the button to record the measurement result into the MEMORY . Press the Button to record the measurement result into the MEMORY ; After selected or the unit will automatically shut-off.
- If you do not press any button, the unit will automatically shut-off after about 2 minutes, and the measurement result will be stored into the MEMORY . To stop the measurement push the START/STOP button.

**Irregular Heartbeat Detector**

Your digital blood pressure monitor features an Irregular Heartbeat Detector. Irregular Heartbeats may influence the results of the measurement. If the monitor detects the Irregular Heartbeats during measurement, the symbol will appear on the display with the measurement values. You can take another measurement to make sure the values are not influenced by moving during measurement or Irregular Heartbeat. If this symbol appears frequently, we recommend you consult your physician. This device does not replace a cardiac examination, but serves to detect pulse irregularities at an early stage.

**Memory Function**

- Recalling Measurements in Memory:
- When the device is off, press the or button to successively view the next previously stored measurements. Measurements are displayed from most recent to the oldest; the memory number appears on the lower left corner.
  - All results for a given measurement are displayed, including measurement results, pulse rate, Hypertension Indicator, and date/time stamp.
  - When the number of readings exceeds 90, the oldest data is replaced with the new record.
  - Press START/STOP button to turn the monitor OFF at any time during review of the stored measurements.

**Clearing Measurements from Memory:**

From power display off, press and hold down the or button until the display shows "CLR" to delete the related memory. This indicates that measurements of or have been erased.

**Pairing up the Blood Pressure Monitor with Your Device Set up process**

- Turn on the smartphone/tablet Bluetooth® connection and Pic Health Station app.
- In the Pic Health Station menu, select "devices" and "add device"
- Press or on your MobileRAPID Wrist.
- The app will automatically detect your MobileRAPID Wrist. Follow the App instruction to complete the pairing.
- During the pairing, only the smartphone/tablet to be connected with the device shall run PIC APP. If you meet problem in pairing please check that no other smartphone in the family are running the APP, or reboot your smartphone.

**Transfer data to your mobile device data transfer during measurement**

After completing the pairing process between the device and the smartphone/tablet just turn on the APP and start the measurement.

After completing the measurement with the APP on, the data are automatically transferred to the active profile.

**Memorized data transfer**

- Turn on the app and select a memorized data in your blood pressure monitor.
- If you press or Each memorized data showed on LCD will be transmitted to your mobile device immediately.

**Maintenance**

- In order to get the best performance from your monitor, please follow the below instructions.
- Put in a dry place and avoid the sunshine.
  - Avoid touching water, clean it with a dry cloth in case.
  - Avoid the intense shaking and collision.
  - Avoid the dusty and unstable temperature environment.
  - Using the wet clothing to remove the dirt.
  - Avoid washing the cuff clean it with a damp cloth

**Calibration and service**

The accuracy of this monitor has been carefully tested and designed for a long life. Indeed the monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or five years of normal use. It is generally recommended to check the unit every two years to ensure correct functioning and accuracy. Please consult the Contact information in this user manual to get in touch with your local authorized Service Center.

- WARNING**
- All the information required for a correct use of the device are included in the instruction for use.
  - The patient is an intended operator.
  - The only operations that are to be considered safe for the patient/user are those listed in the present instruction for use.
  - This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become anesthetic, swollen and even purple due to a lack of blood.
  - This device cannot be used with HF surgical equipment at the same time.
  - This device is not intended for patient transports outside a healthcare facility.
  - This device is intended for adult use only. Keep out of reach of children.
  - People with reduced cognitive abilities shall use the device under supervision.
  - This device is intended for non-invasive measuring and monitoring of arterial blood pressure.
  - It is not intended for use on extremities other than the wrist or for functions other than obtaining a blood pressure measurement.
  - Do not confuse self-monitoring with self-diagnosis. By using this device you can monitor your blood pressure, but it does not replace your doctor. Begin or end medical treatment based solely on medical advice. If you are taking medication, consult your physician to determine the most appropriate time to measure your blood pressure. Never change a prescribed medication without consulting your Physician.

- If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a signal will be displayed. Under this condition, the Electronic Sphygmomanometer can keep function, but the results may not be accurate, it's suggested that you consult with your physician for accurate assessment.
- When the device is used for measurements on patients with common arrhythmias, such as atrial or ventricular premature beats or atrial fibrillation, a deviation may occur. Please consult your doctor about the result.
- If the arm cuff pressure exceeds 300 mmHg it deflates automatically. If the wrist cuff does not deflate when pressure exceeds 300 mmHg, take it off and press the START/STOP button to stop inflation.
- The device is a non-AP/AGP device and not suitable for use in the presence of anesthetic mixtures that are inflammable with air, oxygen or nitrous oxide.
- The operator shall not touch output of batteries and the patient simultaneously.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal.
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- This device is contraindicated for any female who may be suspected of, or is pregnant. Besides provided inaccurate readings, the effects of this device on the fetus are unknown.
- Please use the device under the environment which was provided in the instruction manual. If this is not complied with, correct operation and life of the device may be compromised.
- The materials of the cuff have been tested and found to not cause any potential sensitization or irritation reaction.
- Do not apply the cuff over a wound, on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt is present, on the arm on the side of a mastectomy or lymph node clearance to prevent injury to the patient.
- Do not use more than one ME equipment at the same time on the same limb.
- Please use ACCESSORIES and detachable parts specified/authorized by Manufacturer. Otherwise, it may cause damage to the unit or danger to the user/patients.
- We recommend checking its performance every two years or after it has been repaired. The device must be checked again in case of damage due to impact (such as falls) or exposure to liquids and/or extreme temperatures (hot/cold) or extreme humidity variations.
- This appliance is not intended for public use.
- Please dispose ACCESSORIES, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- Manufacturer will make available to provide technical documentation to Authorities on request.
- In case of problems concerning settings, maintenance or use of the device, contact your local distributor or the Italian toll-free number 800 900 080. Don't open or repair the device by yourself.
- Please use the soft cloth to clean the whole unit. Don't use any abrasive or volatile cleaners.
- The only maintenance allowed is battery change and device cleaning. In order to avoid interruption of the measurement do not perform these activities in use

**Troubleshooting**

PROBLEM	PROBABLE CAUSE	RECOMMENDED ACTION
Nothing appears in the display even when power is turned on	Batteries are drained Batteries are not correctly aligned with terminals	Replace all batteries with new ones Reinsert batteries in the correct position
Low Battery Symbol	Batteries are drained	Replace all batteries with new ones
Device operation time is inconsistent	In colder temperatures batteries have weaker electrical charges Different battery brands have different life spans	Warm up batteries, or use the device in a warmer setting Use Alkaline batteries and replace all batteries at the same time with same brand batteries.
No reading after measurement	Batteries are drained	Replace all batteries with new ones
Suspicious blood pressure results	Perhaps the cuff was improperly positioned Blood pressure naturally varies throughout the day	Adjust patient and wrist cuff to measure Rest a while, relax and measure again
Suspicious heart rate results	Bodily movement during device use Measurement shortly after exercise or exposure to the outdoors	Refrain from moving during measurement. Do not take measurements after exercise or coming back from the outdoors
Power switches off automatically	System design	Push the power button again, and then begin measure again

ERR CODE	MEANING	CORRECTIVE ACTION
ERR 0	No pulse or detect pulses not enough	Take off heavy clothes and retry again.
ERR 1	Leakage in Cuff Pressure/Inflation too low	The wrist cuff is not fastened properly. Re-apply the cuff and take a measurement again.
ERR 2	Measured abnormal values	Take a rest for 3-5 minutes. Re-apply the cuff and take a measurement again.
ERR 3	Inflation fault	Re-apply the cuff and take a measurement again.
	Low batteries	Replace all batteries with new ones.

**TECHNICAL SPECIFICATION**

- Measuring range:** Blood Pressure: 40mmHg ~ 255mmHg  
Pulse Rate: 40 ~ 199beats/min  
Blood Pressure: ± 3mmHg  
Pulse rate: ± 4% of reading  
2 x 1.5V AAA batteries
- Calibration Accuracy:** 90g (exclude batteries)  
91mm x 74.5mm x 15mm (W x H x D)  
approx. 12.5cm ± 21cm
- Power Source:**
- Weight: approx.**
- Dimensions: approx.**
- Cuff circumference (M Size):**
- Operating environment**
- Temperature:** 10°C ~ 40°C (50°F ~ 104°F)
- Relative humidity:** 15% ~ 85% (non condensing)  
700 hPa ~ 1060 hPa
- Atmospheric pressure:**
- Storage/Transportation environment**
- Temperature:** -20°C ~ 50°C (-4°F ~ 122°F)
- Relative humidity:** 15% ~ 85% (non condensing)  
700 hPa ~ 1060 hPa
- Atmospheric pressure:**
- Cuff duration:** 1000 measurement cycles about 250 measurement cycles
- Battery duration:** 5487/8801
- Bluetooth Module N.:** 2.402 ~ 2.435 GHz
- Frequency range:** 3.3 V
- Supply voltage:** 3.3 V
- Maximum output power:** 100mW

**EMC (electromagnetic compatibility) Tables**

Medical Electrical Equipment needs special precautions regarding EMC and needs to be installed and put into service according to the following EMC information. Electromagnetic/Radio Frequency Interference (EMC/RFI): readings may be affected if the unit is operated within a radio frequency electromagnetic field strength exceeding 10 V/m, but the performance of the instrument will not be permanently affected. Care should be taken to keep the blood pressure monitor at least 1 m away from RF transmitters to avoid radio frequency interference.

Guidance and manufacturer's declaration - Electromagnetic emissions		
Phenomenon	Professional healthcare facility environment <sup>a)</sup>	HOME HEALTHCARE ENVIRONMENT <sup>b)</sup>
Conducted and radiated RF EMISSIONS	)	CISPR 11 Group 1 Class B
Radiated RF EMISSIONS (considering the appliance as RF transmitter)	)	ETSI EN 301 489-1 ETSI EN 301 489-17 <sup>c)</sup>
Harmonic Distortion and Voltage Fluctuations and Flicker	Not Applicable	
a) The equipment is suitable for use in Home Health Environment and Professional Health Care Environment limited to patient rooms and outpatient therapy/diagnosis facilities in hospital or clinics. The more restrictive acceptance limits of Group 1 Class B (CISPR 11) have been considered and applied. According to clause 7.1.6 of EN 60601-1-2 the equipment is exempted to meet the emission limit of CISPR 11 in its transmission band [24 ~ 2.4835 GHz]. In the transmission band the emission limits are specified by EN 300 328.		
Guidance and manufacturer's declaration - Electromagnetic immunity - Enclosure port		

Phenomenon	Basic EMC standard or test method	Immunity test levels	
		Professional healthcare facility environment	Home healthcare environment
ELECTROSTATIC DISCHARGE	IEC 61000-4-2	± 8kV contact ± 2 kV, ±4kV, ±8 kV, ±15 kV air	
Radiated RF EM fields	IEC 61000-4-3	)	10 V/m <sup>a)</sup> 80MHz ~ 2.7 GHz 80% AM at 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3		<b>COMPLY</b> NOTE: Further information on separation distances between portable and mobile RF communications equipment (transmitters) and the MobileRAPID Wrist (REF 02 010292 000 000) can be requested to PIKDARE S.r.l., through the contact indicated in the present manual. However, it is advisable to keep the blood pressure monitor at adequate distance of, at least, 1 m from mobile phones or other RF communications transmitters in order to reduce possible interferences.
RATED power frequency magnetic fields:	IEC 61000-4-8	30 A/m <sup>b)</sup> 50 Hz or 60 Hz	
a) The equipment is suitable for use in Home Health Environment and Professional Health Care Environment limited to patient rooms and outpatient therapy/diagnosis facilities in hospital or clinics. The more restrictive IMMUNITY acceptance limits) have been considered and applied. b) Before modulation is applied. c) This test level assumes a minimum distance between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic field of at least 15 cm.			

**Explanation of Symbols**

	"PLEASE REFER TO THE INSTRUCTIONS FOR USE" symbol		"MANUFACTURER" symbol
	Follow instructions for use		Batch code
	Direct current		Catalogue Number
	"TYPE BF APPLIED PART" symbol		Storage and transport temperature symbol
<b>IP22</b>	This product meets the basic safety and essential performance requirements indicated in the IP22 conditioning test (protection against solid foreign objects of 12.5mm Ø and greater and against vertically falling water drops when enclosure tilted up to 15°).		Storage and transport humidity symbol
			Atmospheric pressure
	The CE marking with the Registration Number of the Notified Body. This denotes the compliance of European Medical Device Directive 93/42/EEC.		The Bluetooth® combination mark, Bluetooth® version v. 4.0
	Compliance to RED directive 2014/53/EU		Caution

**Date of manufacture:**

The date of manufacture of the product may be found from the batch coding: 10YYMMXXXX (where: YY=year, MM=month; XXXX=progressive number).

**THE PRODUCT COMPLIES WITH:**

- EN 1060-1:1995/A2:2009 Non-invasive sphygmomanometers – Part 1: General requirements
- EN 1060-3:1997/A2:2009 Non-invasive sphygmomanometers – Part 3: Supplementary requirements for electro-mechanical blood pressure measuring systems
- EN 1060-4:2004 Non-invasive sphygmomanometers - Part 4: Procedures to determine the overall system accuracy of automated non-invasive sphygmomanometers

**THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.**

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

**THIS PRODUCT COMPLIES WITH EU DIRECTIVE 2006/66/EC**

The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for bringing batteries, at the end of their life, to the appropriate collection facilities in order to facilitate treatment and recycling. The correct disposal of your old toy will help prevent potential negative consequences for the environment and human health and favour the recycling of the materials of which the product is made. Illegal disposal of the product by the user will damage the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

**EU DECLARATION OF CONFORMITY:**

Hereby, PIKDARE S.r.l. declares that the radio equipment MobileRAPID Wrist (REF 02 010292 000 000) is in compliance with Directive 2014/53/EU. The full text of the EU declaration of conformity is available at the following internet address: <http://web.pic-solution.com/pic/healthstation/ita/mobile-rapid-wrist/pdf/MobileRAPIDWrist-REF-02010292000000.pdf> In accordance with European Commission Decision No. 2000/299/EC dated 06/04/2000, the frequency band used by this product has been harmonised for all EU countries and this is therefore a Class 1 product which can be used in all European Union countries.

**NOTE:**

- The Bluetooth® trademark and logos are registered trademarks of SIG. Inc. and any use thereof by PIKDARE is under license.
- Google Play and the Google Play logo are Google Inc. trademarks.
- Apple and the Apple logo are trademarks of Apple Inc., registered in the United States and in other countries. App Store is a service trademark of Apple Inc. registered in the United States and in other countries.

**Pouch cleaning symbols**

- Wash by hand in cold water,
- Do not bleach,
- Do not tumble dry,
- line drying,
- Do not iron,
- Do not dry clean.

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