










Mister 8

CE
0068

ISTRUZIONI D'USO

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www.picsolution.com

INSTRUCTIONS FOR USE

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer,

Thank you for choosing PIC Solution's Mister 8, a micro piston aerosol therapy unit that was designed and manufactured using state-of-the-art technology.

The Mister 8 is a class IIa medical device that converts medicines from liquid or suspension form to aerosol form, so that they can be administered directly to the airways.

The Acti-Fast Pro nebuliser cup allows for deeper penetration of the medicine into the respiratory tract, thanks to the optimal size of the aerosol particles.

The Mister 8 was designed to contain all the accessories that come with it (nebuliser cup, mouthpiece, nose fork, mask for adults and children).

Before using the device, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

For your convenience, consult the Quick User Guide found inside the package.

ACCESSORIES

A - Adult mask

B - Children's mask

C - Mouthpiece

D - Nose fork

E - Air tube

F - Pic Acti-Fast Pro nebuliser cup

F1 - Aerosol cup junction - accessories

A, B, C, D, E, F, F1, I, are available as replacement accessories in the AIRFamily evolution kit  02038405200000

DEVICE COMPONENTS

G - On/Off button

H - Compressed air outlet

I - Filters

L - Nebuliser cup support

M - Octopus accessories storage container

N - Mains adapter connector

GENERAL WARNINGS

- ⚠️ **WARNING:** In spite of its shape, this device is not a toy. This product must be assembled, maintained and handled by an adult. Never leave children alone during use: the device contains small parts that could be swallowed.
- This product must remain in place for use only for the time strictly needed. Once treatment has been completed, the device must be stored in a safe place out of the reach of children.
- Before use, carefully read the information in this manual, and keep for future reference. If you have any doubts about the meaning of the contents of this instruction booklet contact the dealer, distributor or manufacturer.
- Device designed for aerosol therapy. Always follow the instructions of a medical practitioner as to the type of medicine to use, the dose, and the frequency and duration of inhalations. Any use other than the one for which the device is intended is considered improper and therefore dangerous; PIKDARE S.r.l. cannot be held liable for any damage caused by incorrect and/or unreasonable use or if the device is used in electrical systems that do not conform to applicable safety regulations.
- The device may stop working in the event of a power outage, or following sudden breakdown or adverse conditions of any kind; you are therefore advised to have an alternative medicine or device on hand (agreed upon with your medical practitioner) for use in such situations.
- After removing the device from its packaging, check that it is intact, without any visible damage that may have occurred during transport. If in doubt, do not use the device, and contact the PIC Service Centre (www.picsolution.com).

- Do not plug in or unplug the AC/DC adapter and/or touch the pin connector (metal part) of the AC/DC adapter with wet hands.
- IMPORTANT: to connect the device to the mains, please connect the supply pin connector to the device supply socket (N) and only then connect the AC/DC adapter plug to the mains socket. To disconnect the device from the mains, first disconnect the AC/DC adapter plug from the mains socket and only after that disconnect the pin connector from the device supply socket (N).
- Packaging (bags, box, etc.) may be dangerous and must be kept out of the reach of children. Keep the device out of reach of children, pets and persons with cognitive limitations.
- Before connecting the mains adaptor to the mains socket, ensure that the information on the power data plate and the mains adaptor are compatible with your electricity distribution network and the corresponding socket. The reference power data plate can be found on the mains adaptor.
- Should the AC/DC adapter plug supplied be incompatible with the electrical socket, use certified plug adapters, provided applicable legislation in the individual country allows for it.
- To prevent overheating and damage to the compressor, switch off the device for at least 40 minutes after every 20 minutes of uninterrupted use.
- ⚠️WARNING! Using the device before 40 minutes have elapsed may result in the motor overheating and may cause the safety cut-off device to activate.
- To ensure correct operation, make sure that the filter is dry. The filter cover may be removed using a tool with a flat blade. ⚠️ WARNING! Keep the filter cover and the filter out of the reach of children as these are small parts and could represent a swallowing and suffocation hazard.
- Do not use the device around nitrous oxide, oxygen or anaesthetic mixtures that are flammable in air.
- Turn off the device and disconnect the AC/DC adapter from the mains after use and/or before adding more medicine. Do not fill the nebuliser cup past the maximum level indicated (10 ml) and check that all components have been assembled correctly.
- Keep the device, the AC/DC adapter and the power cord away from hot surfaces.
- Do not use the device while in the bath or shower, or in damp environments, or near baths, sinks, washbasins, or in any other situation in which there are liquids that may come into contact with the device.
- Never touch the device and the AC/DC adapter with wet or damp hands.
- Do not drop or put the device into water or other liquids. Should the device fall into water, immediately unplug the device, do not use the device anymore and contact the PIC Service Centre (www.picsolution.com).
- Do not use the device if you feel drowsy.
- If children or persons with physical or mental disabilities use the device, proper supervision by an adult is required. Children must be kept under adult supervision to ensure that they do not play with the device. The device contains small parts that can be swallowed and the AC/DC adapter power cord is a strangling hazard.
- Do not use connectors or accessories not recommended by the manufacturer.
- Make sure the device is placed on a flat, level surface when in use in order to prevent spillage.
- During use, the surface must be free of any objects that may obstruct the correct flow of air.
- Do not obstruct the holes of the filter cover, which are located in the bottom of the device, with cloths or anything else; the motor may overheat, reducing the performance of the device.
- Do not pour liquids on the device, especially on the filter cover, since this may lead to the undesired expulsion/nebulisation of those liquids.
- Before cleaning or performing any maintenance work on the device, disconnect the AC/DC adapter from the mains by pulling the plug out of the electrical socket.
- If you decide you will no longer use the device, disconnect the AC/DC adapter from the mains and cut off the power cord so that the device cannot work. Dispose of the device and the AC/DC adapter in compliance with applicable laws. We also recommend eliminating any parts that are potential hazards, particularly for children.
- Refer to applicable laws for information about the disposal of accessories that are subject to wear. To dispose of the device, refer to Directive 2012/19/EU.
- Do not use the device and the AC/DC adapter if, after a fall, there are signs of damage or deterioration to any of its components. If in doubt, do not use the device or the AC/DC adapter, and contact the PIC Service Centre (www.picsolution.com).
- In the event of a breakdown and/or anomaly when using the device, immediately turn it off and disconnect the AC/DC adapter plug from the mains socket. Then contact the PIC Service Centre (www.picsolution.com) exclusively.
- ⚠️WARNING! The AC/DC adapter is the element that separates the device from the electrical mains. When using the device, to disconnect it from the mains, switch off the device and remove the AC/DC adapter plug from the power outlet. Only connect the AC/DC adapter plug to a power socket that is easy to reach, so that it will be easy to disconnect the device from the mains. Do not attempt to open, repair or make changes to the device and/or the AC/DC adapter.
- The patient and operator are one and the same individual. The only operations that are to be considered safe for the patient/user are those listed in the following paragraphs: "Preparing and using the device", "Cleaning and maintenance of the compressor", "Cleaning and maintenance of the nebuliser cup and accessories".

- ⚠️**WARNING!** Do not service and/or repair the appliance while it is being used! Any other action (including disassembling and later reassembling the device) or any maintenance/service operations must be done solely by a PiC Service Centre.

PREPARING AND USING THE DEVICE

- Take the device and the accessories out of the package and check that they are intact; in the event that they have been damaged, replace them or refer to a PiC Service Centre (www.picsolution.com);
- Add the medicine and/or the saline solution to the nebuliser cup reservoir, using the dose recommended by your physician (MIN 2ml - MAX 10ml);
- Close the nebuliser cup;
- Add the nebuliser cup junction (F1) to the upper cylinder of the nebuliser cup;
- Connect the air tube to the nebuliser cup;
- Add the accessory you plan to use for the treatment: mask for adults or children, nose fork, or mouthpiece (Applied parts made with biocompatible material) to the nebuliser cup junction. In case of doubt, ask your physician about which accessory to use.
- Connect the tube to the compressed air outlet (H);
- Connect the AC/DC adapter cable pin connector to the device socket (N).
- Connect the AC/DC adapter plug to the mains power socket;
- Switch ON by pressing the ON/OFF button and start the treatment;
- Perform treatment preferably while seated in a comfortable position; when the aerosol jet becomes intermittent, interrupt the treatment for a few seconds and let the suspended drops of solution settle on the sides of the nebuliser cup. Resume the treatment session, which will be finished when no more nebulised solution is released from the nebuliser cup;
- Once the treatment session is completed, switch the device OFF, unplug the AC/DC adapter plug from the mains socket, unplug the AC/DC adapter cable pin connector from the device socket (N), remove the accessory used and clean the device and accessories according to the instructions in the "Cleaning and Maintenance" paragraph.

CLEANING AND MAINTENANCE

COMPRESSOR CLEANING AND MAINTENANCE

⚠️**WARNING!** Do not wet the device or put it into water or other liquids. For cleaning, use a clean, dry cloth only.

CLEANING AND MAINTENANCE OF THE NEBULISER CUP AND ACCESSORIES

Remove and detach all the nebuliser cup components and submerge them in boiling water for about 5 minutes.

⚠️**WARNING!** Do not boil the masks.

- Rinse the accessories (mask, mouthpiece, nose fork and tube) under running water.
 - Cold sterilise all parts that come in contact with the patient using denatured alcohol.
 - Carefully dry all accessories and the device before storage.
 - Store the device and all accessories in a clean, cool, dry place, away from light and heat.
 - Never use benzene, thinners or other flammable chemical substances for cleaning.
 - For hygienic purposes, never use the same accessories on more than one patient. Purchase a dedicated kit for each user (AIRFamily evolution kit [\[REF\] 02038405200000](#)).
 - Clean and check the condition of the accessories before and after use. Replace them if damaged.
- ⚠️**WARNING:** Do not leave accessories unattended since small parts of the nebuliser cup, mouthpiece and nose fork may be ingested by children and are therefore a choking and suffocation hazard. The air tube may be a strangling hazard.

CHECKING AND REPLACING THE FILTER

Periodically check the condition of the filter. The filter was inserted to protect the compressor.

Correct filter maintenance prolongs the life of the aerosol. The frequency with which filter should be replaced depends on the conditions in which the device is used.

The filter is located on the bottom of the aerosol.

TO REPLACE THE FILTER:

- Open the filter cover cap with the help of a tool;
- Remove the filter to be replaced from its housing
- Carefully fit the new filter

- Fit the filter cover cap back on.

⚠WARNING! Never leave the filter cover cap and/or filters unattended: These are small parts posing the risk of suffocation if swallowed by children.

TROUBLESHOOTING

PROBLEM	CAUSE	SOLUTION
Scarce nebulisation	Nebuliser cup clogged	Remove and disconnect all nebuliser cup components and clean them according to the instructions provided in this instruction manual.
Scarce nebulisation	Nebuliser cup clogged	If the problem remains after cleaning, replace the nebuliser cup.
No nebulisation	Components of the nebuliser cup assembled incorrectly	Disassemble and reassemble the nebuliser cup. If the problem remains, replace the nebuliser cup.
No nebulisation	The compressor does not work	Check whether the AC/DC adapter plug is disconnected from the mains socket, check that the ON/OFF button is ON, that the mains socket is operational, and that the AC/DC adapter power supply cable pin is connected to the device socket (N).
Slow nebulisation	Medicine too thick	Dilute the medicine as specified by the pharmaceutical company.

TECHNICAL SPECIFICATIONS

External power supply:	100-240 VAC 50/60Hz 0.6A 12V 2A REF 02010298000000 (EU plug) REF 02010299000000 (UK plug)
Service life in storage conditions:	5 years
Service life of accessories:	1 year
Frequency:	50/60 Hz
Maximum flow:	12 ± 2l/min
Operational flow:	5 ± 1l/min
Max pressure:	2 ± 0.4 bar
Operational pressure:	0.6 ± 0.18 bar
Conditions of use:	20' ON – 40' OFF at 40°C
NEB RATE:	> 0.30 ml/min
MMAD:	< 3 µm
Weight of the Mister 8 [without accessories]:	~800g
Nebuliser cup capacity:	Max 10ml; Min 2ml
Breathable fraction:	> 75%

Device service life when in use: 1000 cycles/treatments

Noise level at 50cm: < 56,5dB(A)

The information provided by the manufacturer on the device's performance capabilities in accordance with EN 13544-1 may not apply to medicines provided in suspension form or those that are highly viscous.

EN 60601-1 Medical electrical equipment – Part 1: General requirements for safety.

EN 60601-1-2 Medical electrical equipment – Part 1: General requirements for basic safety and essential performance.

- Collateral standard: electromagnetic compatibility - Requirements and tests.

EN 13544-1 Respiratory therapy equipment – Part one: Nebulising systems and their components

This section contains information specific to product compliance with the EN 60601-1-2 standard. The Mister 8 is a medical electrical device that requires special precautions regarding electromagnetic compatibility and needs to be installed and commissioned in accordance with the electromagnetic information provided.

Mobile and portable RF communications equipment (mobile phones, transceivers, etc.) may affect the device.

The Mister 8 [REF] 02038211000000 was designed to be used in an electromagnetic environment in accordance with the specifications provided below. The person who buys or uses the Mister 8 [REF] 02038211000000 must ensure that the device is used in an environment that complies with these specifications.

Guidelines and Manufacturer's Declaration – Electromagnetic Emissions

Phenomenon	Professional healthcare facilities ^{a)}	HOUSEHOLD ENVIRONMENT ^{a)}
Conducted and radiated radiofrequency EMISSIONS	^{a)}	CISPR 11 Group 1 Class B
Harmonic distortion	IEC 61000-3-2 ^{b)} Class A	
Voltage fluctuations/flicker	COMPLIANT WITH IEC 61000-3-3 ^{b)}	

^{a)} The device is destined for household use or use in professional healthcare facilities and may only be used in the patient's room and in locations designated for respiratory therapy in a hospital or clinic. The most restrictive acceptance limits established for Group 1 Class B (CISPR 11) have been applied. The device is suited for use in the aforementioned environments so long as it is connected to the public mains.

^{b)} The test may be applied to this environment provided that the ME EQUIPMENT and the ME SYSTEM used are connected to the PUBLIC MAINS and that the power supply is compliant with the prescriptions set by the reference policy governing electromagnetic compatibility (EMC).

Guidelines and Manufacturer's Declaration - Electromagnetic Immunity - Enclosure Door

Phenomenon	EMC Regulation or test method	Immunity test levels	
		Professional healthcare facilities	Household environment
ELECTROSTATIC DISCHARGES	IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	
Radiated, radio-frequency, electromagnetic fields	IEC 61000-4-3	^{a)}	10 V/m ^{b)} 80 MHz - 2.7 GHz 80% AM at 1 kHz

Fields located near RF wireless communication devices	IEC 61000-4-3	COMPLIANT NOTE: Further information about distances to be maintained between portable and mobile RF communications equipment (transmitters) and the Mister 8 REF 02038211000000 can be requested from PIKDARE S.r.l. using the contact information provided in this manual. It is nevertheless advisable to keep a safe distance (at least 1 m) between the electromechanical aerosol therapy device and mobile phones or other RF communications devices in order to minimise possible interference.
Magnetic fields at NOMINAL mains frequency.	IEC 61000-4-8	30 A/m ³) 50 Hz or 60 Hz

^a) The device is destined for household use or use in professional healthcare facilities and may only be used in the patient's room and in locations designated for respiratory therapy in a hospital or clinic. The most restrictive limits were applied for the immunity tests.

^b) Before modulation.

^c) This test level presupposes a minimum distance of at least 15 cm between the piece of piece of ME EQUIPMENT or the ME SYSTEM and the sources of electromagnetic fields at mains frequency.

Guidelines and Manufacturer's Declaration - Electromagnetic Immunity - Alternating Current Entry Port

Phenomenon	EMC Regulation or test method	Immunity test levels	
		Professional healthcare facilities	Household environment
Electrical fast transient/burst	IEC 61000-4-4	± 2 kV 100 kHz frequency of repetition	
Line-to-earth surges	IEC 61000-4-5	± 0.5 kV, ±1 kV	
Line-to-earth surges	IEC 61000-4-5	Class II device or equipped with Class II AC/DC adapter. Test not applicable	
Conducted disturbances induced by radio frequency fields.	IEC 61000-4-6	3 V ^a) 0.15 MHz - 80 MHz 6 V ^a) for ISM bands between 0.15 MHz and 80 MHz ^b) 80% AM at 1 kHz	
Voltage dips.	IEC 61000-4-11	0% in UT; 0.5 cycles At 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°.	
		0% in UT; 1 cycle and 70% in UT; 25/30 cycles Monophase: at 0°	
Short interruptions ^c)	IEC 61000-4-11	0% in UT; 250/300 cycles	

^a) RMS value before modulation.

^b) The frequency bands from 0.15 MHz to 80 MHz used by industrial, scientific and medical (ISM) radiofrequency devices are: 6.765 MHz - 6.795 MHz, 13.553 MHz - 13.567 MHz, 26.957 MHz - 27.283 MHz and 40.66 MHz - 40.70 MHz. The radio bands from 0.15 MHz to 80 MHz used for amateur radio are: 1.8 MHz - 2.0 MHz, 3.5 MHz - 4.0 MHz, 5.3 MHz - 5.4 MHz, 7 MHz - 7.3 MHz, 10.1 MHz - 10.15 MHz, 14 MHz - 14.2 MHz, 18.07 MHz - 18.17 MHz, 21.0 MHz - 21.4 MHz, 24.89 MHz - 24.99 MHz, 28.0 MHz - 29.7 MHz and 50.0 MHz - 54.0 MHz.

^c) ME EQUIPMENT with direct current input for use with an AC/DC adapter, tested using an adapter compliant with the specifications provided by the ME EQUIPMENT MANUFACTURER. THE IMMUNITY TEST LEVEL was applied at the adapter's alternating current entry point.

LEGEND OF SYMBOLS:

Warning! Follow the Instructions for Use

I-O

ON/OFF

Double insulation

Hz

Frequency

IP21

Protected against the ingress of solids and liquids. Protected against dripping water.

SN

Serial Number

Device with type BF components

REF

Product ID code

Complies with Directive MDD 93/42/EEC and 2007/47/EC

Manufacturer

Alternating current

LOT

Batch

Direct current

DEVICE DATE OF PRODUCTION: the second pair of digits in the lot number identifies the year, the third pair identifies the month. (**LOT** 1016100001; 16 = 2016; 10 = October)

*Data on noise levels were measured on a new device. Values may vary over time with use.

CONDITIONS FOR STORAGE AND TRANSPORTATION:

Relative humidity: 15% / 85% RH

(store in a cool, dry place)

Temperature: -25°C / 70°C

Atmospheric pressure: from 700hPa to 1060hPa

ENVIRONMENTAL CONDITIONS DURING USE:

Relative humidity: 15% / 85% RH

Temperature: 5°C / 40°C

Atmospheric pressure: from 700hPA to 1060hPA

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your retailer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the unwanted appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

WARRANTY

The product is protected by a 5-year warranty against any conformity defect in normal conditions of use as provided by the instructions for use.

Therefore, the warranty will not apply in the event of damage caused by improper use, wear or accidental events. For the duration of the warranty against conformity defects please refer to the specific provisions of national laws.

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