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Rev. 02-10/2023
76010394-000-000_02_10_2023
 02010394000000

THESE INSTRUCTIONS ARE IMPORTANT. PLEASE KEEP THEM FOR FUTURE REFERENCE.

Dear customer,
Thank you for choosing Pic Solutions OneRAPID, an automatic blood pressure meter designed and manufactured according to the most up-to-date technologies. Before using the device, we highly recommend that you read the short list of warnings contained in this manual, to make sure you have correctly understood how to use it.

INTRODUCTION

The regular use of a Blood pressure monitor allows you to monitor blood pressure daily, but it should not take the place of regular check-ups by your doctor. We recommend that you visit your doctor regularly for general check-ups and for more detailed information on blood pressure. Blood pressure measurements determined with this unit are equivalent to those obtained by a trained observer using cuff/stethoscope auscultation method, within the limits prescribed by the international (EN ISO 81060-2:2019+A1:2020) validation protocol.

WARRANTY

The product is guaranteed for 5 years against all conformity defects under the normal utilization conditions provided in the user instructions. Therefore, the warranty shall not apply in case of damages caused by improper use, wear, or accidental events. Batteries, cuff and wearing parts are not included in the warranty. In case of malfunction, contact exclusively the authorized Pic Technical Customer Service.

INTENDED USE

This device is intended for non-invasive measuring of arterial blood pressure in adults in a home or professional setting. An algorithm is used to derive the heart rate from the blood pressure measurements and to signal the presence of an irregular heart rate.

INTENDED POPULATION

This device is not suitable for use on newborns or infant (up to 3 years age); the device parameters are set for adults only. The device has not been validated for use on pregnant or pre-ecliptic women and may be dangerous. Consult with your physician before using this monitor if you suffer from one of more of the following conditions: common arrhythmias such as premature ventricular beats or atrial fibrillation, peripheral arterial disease (in cases of serious arteriosclerosis - hardening of the arteries), poor perfusion, diabetes, pregnancy, pre-eclampsia, renal diseases, trembling, shivering or if you have active implanted medical devices. Please note that any of these conditions may affect measurement readings. Suitability for use on children and adolescents should be assessed by a physician. The device is not suitable for use on persons with arm circumferences outside the range of 22 to 42 cm. Persons with reduced cognitive capacity should use the device under supervision.

RAPID TECH Comfort Experience™

This unit uses the oscillometric method to detect your blood pressure. RAPID TECH Comfort Experience™ is the latest generation technology chosen by Pic for the blood pressure meters of the line RAPID. It allows fast and, therefore, less annoying monitoring because the measurement is performed already during the cuff inflation phase. Pic OneRAPID has been clinically validated in accordance with the international (EN ISO 81060-2:2019+A1:2020) validation protocol, as a guarantee of its accuracy and precision.

NAME/FUNCTION OF EACH PART

ABOUT BLOOD PRESSURE

What is blood pressure?
Blood pressure is the pressure exerted on the artery tube while blood flows through the arteries. The pressure measured when the heart contracts and sends blood out of the heart is systolic (highest) blood pressure. The pressure measured when the heart dilates with blood flowing back into the heart is called diastolic (lowest) blood pressure.

Why measure your blood pressure?


Among the various health problems afflicting modern people, problems associated with high blood pressure are by far the most common. High blood pressure's dangerously strong correlation with cardiovascular diseases and high morbidity has made measuring blood pressures a necessity of identifying those at risk.

Blood Pressure Fluctuation

Blood pressure fluctuates all the time! You should not be over worried if you encountered two or three measurements at high levels. Blood pressure changes over the month and even throughout the day. It is also influenced by the season, temperature and activities.

Classification of high blood pressure

The out-of-office blood pressure classification published by ESH/ESC in 2018 is as follows:

Note: It is not intended to provide a basis of any type of rush toward emergency conditions/diagnosis based on the scheme: the scheme is meant only to discriminate between the different levels of blood pressure. Light on START/STOP button shows the blood pressure level according to ESH classification: SYS<135mmHg and DIA<85mmHg, the light is green - SYS>=135mmHg or/and DIA>=85mmHg, the light is orange.

WARNING

Only a physician can tell your normal BP range. Please contact a physician if the average of your measuring results falls out of the range. Kindly note that only a physician could tell whether your blood pressure value has reached a dangerous point.

Irregular Heartbeat Detector: This blood pressure meter is equipped with an algorithm for Irregular Heartbeat (IHB) Detection. An irregular heartbeat is detected when a heartbeat rhythm varies while the unit is measuring the systolic and diastolic blood pressure. During each measurement, the monitor records all the pulse intervals and calculate the average; if there are two or more pulse intervals, the difference between each interval and the average is more than the average value of ±25%, or there are four or more pulse intervals, the difference between each interval and the average is more than the average value of ±15%, then the irregular heartbeat symbol will appear on the display with the measurement result.

WARNING

The appearance of the IHB icon indicates that a pulse irregularity consistent with an irregular heartbeat was detected during measurement. Usually this is NOT a cause for concern. However, if the symbol appears often, we recommend you seek medical advice. We wish to underline that the device does not replace a visit to the cardiologist, but it is useful to detect any irregularities in early stage.

Why the blood pressure I get from the hospital is different from home?

The blood pressure is different even during 24 hours because of the weather, emotion, exercise etc., specially the "white coat" in hospital which makes the results are higher than the ones at home.

If the result is the same if measuring on the right arm?

Measurements can be taken on either arm; results can however differ sometimes, so it is always advisable to use the same arm.

Display explanation

At starting check the completeness of the segments

SYMBOL	DESCRIPTION	EXPLANATION
188	Memory group number	Share the segments with heart rate "188" number.
[Battery icon]	Low power symbol	Running out of power and needs to be recharged asap.
[Person icon]	Movement symbol	Indicates body movement was detected.
OK	Cuff tightness detection symbol	When only the 1st icon show up, cuff tightness detection on-going. When both icons light-up, detection is done and cuff tightness OK.
[Bluetooth icon]	Bluetooth® symbol	When it stop blinking, it means Bluetooth® connected.
[Book icon]	Book symbol	Memory data need to be uploaded to app.
AVG	Average symbol	Average value of the last 3 measurements within 30min.
1 2	User symbol	User 1, user 2
mmHg	mmHg	Blood pressure measurement unit, mmHg
[Heart icon]	Heartbeat indication	Heartbeat signals have been detected.
[IHB icon]	IHB symbol	Irregular heartbeat was detected.
SYS	SYS	Systolic pressure
DIA	DIA	Diastolic pressure
PULSE /min	PULSE	Heart rate per minute
MEM	Memory symbol	Indicates entered into memory query mode.

Power supply

The battery of the OneRAPID is a built-in rechargeable Li-polymer battery type, the battery is 3,6 V 1000 mAh.

Charging the power under following circumstances:

- Is displayed the warning voltage symbol
- The LCD display dims,
- When powering on the monitor, the display does not light up.

WARNING VOLTAGE	LOW VOLTAGE
[Warning Voltage Symbol]	[Low Voltage Symbol]
If the battery is low, the low battery symbol is displayed. This warning has the highest priority. For example, the low battery symbol should be displayed in set mode, query mode, pairing mode, etc., or other error code states	If the battery is extremely low, the low battery symbol and "bAt Lo" will be displayed, and it will shut down automatically after 5s

On any interface, when the device is at the warning voltage, the symbol of low power will be displayed continuously until the device is charged with AC/DC adaptor. Please use only the AC/DC adaptor (REF) 02010309100000 (EU) or (BLJ06L050100U-V) (REF) 02010308100000 (BLJ06L050100U-B) (UK) to charge the battery. If there is overvoltage, the prompt code "USB Er" will continue to display on the screen. After 5s, the device will automatically shut down.

In any display interface, when entering the warning voltage condition, the low power symbol continues to display, while the content under the current interface remains unchanged.

NOTE: while charging the device cannot be used and all the device functions are deactivated. The charging cycle will last up to 3,5 hours.

While charging the colored lighted line between on start/stop button will be blue and steady on; when fully charged, the light goes off. If the battery is fully charged, it can be used for about 100 times. The battery of OneRAPID is built-in rechargeable lithium-ion battery, please do not disassemble it by the unauthorized maintenance personnel. If the battery cannot charge the power normally or the blood pressure monitor cannot use normally, please connect with the authorized maintenance personnel.

Before the first usage - Setting

To turn on the device the first time long press the START/STOP button and it will automatically enter the setting mode.

Specific operation of setting mode

Bluetooth® pairing: After entering the setting mode, access the Pic Health Station APP and be sure that the Bluetooth® of the smartphone is on. The Bluetooth® function will automatically turn on in background and the Bluetooth® symbol will blink. **If the pairing is successful, the device time will be automatically synchronized to the time of the mobile phone end. The device end does not display the time, but has an internal real-time clock running, so as to mark the time stamp for each measurement data.** If a large amount of historical data exists after the device obtains the synchronization time from the mobile phone, it needs to automatically recall the measured historical data based on the time difference and reassign the time stamps to each group of historical data. After the completion of Bluetooth® pairing, the Bluetooth® symbol will be displayed continuously, and the device will automatically shut down after the screen displays "done". If the pairing procedure fails the message "Er" will be displayed and the device will shut down. If you want to interrupt the setting at any step during the setting process, hold down the START/STOP button or press no button for 30 seconds before the system automatically shuts down.

Bluetooth® Chip No. : BK3432
Frequency Range: 2402Mhz-2480Mhz
Supply Voltage: 3.3V
Maximum Output Power: 4dB
Transmitting distance: 10 m

User change

In the shutdown state, press the START/STOP button on the device to start measurement, and at the same time, you can perform user switch: press User selection button to change the user. After you press the START/STOP button to power on, 1 or 2 are displayed. By default, the last user blinks, indicating that you are under the current user. Press the user selection button once more to switch the user before the measurement starts. If you wait for more than 2 seconds, the current user will be automatically confirmed and the measurement will begin immediately. In the shutdown state, press the users' selection key to switch the user.

Helpful Tips

- Here are a few helpful tips to help you obtain more accurate readings:
- Do not measure your blood pressure immediately after consuming a large meal.
 - To obtain more accurate readings, please wait one hour after meal before measuring.
 - Do not smoke or drink tea, coffee or alcohol before measuring your blood pressure.
 - Do not measure your blood pressure immediately after taking a bath: wait at least 20 minutes.
 - You should not be physically tired or exhausted while taking a measurement.
 - Do not measure your blood pressure when you want to discharge urine.
 - It is important that you relax during measurement.
 - Do not take measurements if you are under stress or tense.
 - Take your blood pressure at normal body temperature. If you are feeling cold or hot, wait a while before taking a measurement.
 - Rest for 5 minutes before measuring.
 - Carry out 3 measurements in 3 minutes. Do not exceed the three sequential measurements to avoid side effect on normal blood circulation.
 - Take the measurement in a silent room.
 - Do not move and do not talk when taking a measurement
 - Do not cross your legs and keep your feet on the ground.
 - Keep your back against the backrest of the chair.
 - For a meaningful comparison, try to measure under similar conditions. For example, take daily measurements at approximately the same time, on the same arm, or as directed by a physician.

Applying the Cuff

1. Remove all watches, jewellery, etc. prior to attaching the arm monitor. Clothing sleeves should be rolled up and the cuff should be wrapped on bare skin for correct measurements.
2. Sit comfortably and apply the cuff preferably to left arm with palm facing up as Fig. A.
3. Make sure the edge of the cuff is about 2-3 cm from the elbow.
4. In order to ensure accurate measurements, fasten the velcro strap securely around your arm so there is no extra space between the cuff and your skin. If the cuff is too loose, the measurement will not be accurate.
5. Position the device body over the main artery (on the inside of your arm). Note: Locate the main artery by pressing with 2 fingers approximately 2 cm above the bend of your elbow on the inside of your left arm. Identify where the pulse can be felt the strongest. This is your main artery.

6. The cuff should maintain at the same level as the right atrium of the heart.
7. If you apply the cuff to the left arm, as further support to the correct positioning of the device you can use the red band as a reference; apply it over the main artery as described in point 5.

After the cuff detection symbol shows OK, the pressure starts.

Single measurement

In the power off state, press the [Power button/User] to enter the single measurement mode. Short Press "START/STOP" to enter the blood pressure measurement Mode in power off State. All the segments will be showed to verify the completeness of the display.


After you press the START/STOP button to power on, 1 or 2 are displayed. By default, the last Area used blinks, indicating that you are under that user. Press the users' selection button to switch the memory area. If you wait for more than 2 seconds, the current user will be automatically confirmed and the measurement will begin immediately. The measurement starts, please remain still and quiet. The cuff starts to inflate. The APP pairing can be done also during the whole measurement process; in case of APP pairing without data transmission, the Bluetooth® symbol flashes, and the Bluetooth® symbol continues to display after the pairing is successful. Also during the measurement result display phase, Bluetooth® pairing and data transmission can be performed. After successful data transmission, the device will automatically shut down after no operation for 30 seconds. The cuff detection symbol (flashing) and the pressure value are displayed during the inflation process; after the cuff detection symbol shows OK, the inflation starts. When the heart rate is detected, the heartbeat indicator starts to flash, and the flashing frequency is consistent with the pulse. If the cuff detection symbol is not displayed OK, no heartbeat indicator will appear.

If the arm shakes during the measurement, the Irregular Body movement symbol is displayed.

[Cuff icon]	[Cuff icon]	[Pulse icon]
Inflating-cuff tightness detection in progress	The cuff tightness is ok	Pulse detected
[Irregular Body icon]	[Irregular HeartBeat icon]	color indication of the pressure level
Irregular Body movement symbol	Irregular HeartBeat	

If the measurement result is out of the measurement range (SYS: 60 ~ 230mmHg or DIA: 40 ~ 130 mmHg or PULS 40 ~ 199 Pulse/min), the device will display "OUT".

After the blood pressure measurement is completed, the screen will display the results in this way: the value of the systolic and diastolic blood pressure alternately with that of the heart rate, blood pressure level classification light, irregular heartbeat (if any), irregular body movement (if any). The result is stored to the current user and the book symbol will be uploaded is displayed. If the Bluetooth® is connected, the Bluetooth® transmission will start. During the process of Bluetooth® transmission, the screen will be steady on with the Bluetooth® symbol. When the transmission is completed, the book symbol will disappear immediately, and the Bluetooth® symbol will disappear automatically after keeping it steady for 5s. After the measurement result is displayed, press the "Start/stop button" to shut down the device or wait for data upload. After the data upload succeeds or fails, the system automatically shuts down after 30 seconds without any operation. During measurement, if irregular body movement is detected, the screen will blink symbol for 5s. After 5s, the screen will check whether the irregular body movement is detected again. If the arm movement is not detected anymore, the IBM symbol will not be displayed but it will be shown together with the result and stored in the record; if the irregular body movement is so strong that it could affect the measurement accuracy, E2 will be displayed. When data upload fails within 30s, book symbol shows, the Bluetooth® symbol by normally on become a blinking, if in the 30s to establish the connection is successful, the Bluetooth® symbol back to normally on continue to complete the data transmission, if not to establish a connection in the 30 s, are judged to be connection timeout, directly to turn it off. During measurement, when the device detects an error, the screen continues to display the corresponding error code. During measurement, you can press the "START/STOP button" at any time to shut down the device and deflate the cuff at the same time.

Memory query

In the shutdown state, press Memory Button to enter the query mode of user data memory. The user symbol 1 or 2 flashes, the memory symbol "MEM" is continuously displayed. Press the User button to switch the user, if needed. Press Memory button to select and enter directly in the memory query selected. If you wait for more than 2 seconds, the user is automatically confirmed and you enter directly in the selected memory area. Display the average value of the last three measurements, if taken within 30 minutes, with the average symbol AVG. If there are less than 3 records or they are not taken within 30 minutes, the latest measurement record (00 group) will be directly displayed after entering the memory mode. Press the Up key or the Down key to start the query. If the user does not press any key, the user automatically confirms the query after three seconds and starts the query from the first group. For every record, it is first shown the systolic and diastolic blood pressure value and then the heart rate and record group. Memory query display sequence: AVG → 0 group → 1 group → ... → Oldest (199 group) Press the Down key to Increase the number of memory groups (AVG → 0 → 1 → ...); Press

the Up key to reduce the number of memory groups (AVG → oldest in history → ...). The most recent recording (0) is shown first. Every time a new recording is made, the first (0) record is assigned to it. All the others are moved by one figure (e.g. 0 becomes 1, and so on), whereas the last record (199) is deleted from the list. If the user has no record and all the values are displayed as "–"; The memory query process can be shut down at any time by pressing the START/STOP, or automatically shut down after 30s operation without any key.

Delete a single record

After entering the memory query mode, press the Up or Down keys to go to the record you want to delete, and hold down the Down key for 3s. "dEL y" will flash on the screen. You can press the Up button or the Memory button/Down button to switch to "dEL y" or "dEL no". Press the START/STOP button to confirm. If "dEL y" is confirmed, this record is deleted and done is displayed continuously. After 1 seconds, the record is automatically switched to the previous memory data. If "dEL no" is confirmed, the deletion action is cancelled, and the record is retained on the original query page. The single delete mode cannot be entered if there isn't any record stored.

Delete all records of the current user

Enter the memory query mode, hold down the Up and Down at any record for 3s at the same time, and "dEL AL" will flash on the screen. You can press the Up or Down to switch to "dEL AL" or "dEL no". Press the START/STOP button to confirm. If "dEL AL" is confirmed, all records of the current user are deleted and "done" is displayed continuously. After 3 seconds, all values are displayed as "–". If "dEL no" is confirmed, the deletion action is cancelled, and the record is retained on the original query page. The delete all records mode cannot be entered if there isn't any record stored.

ERROR code

SYMPTOM	CAUSE	MEASURE
The display screen is not lit.	The battery is exhausted.	Need charging
	The power adapter plug is not plugged in.	Plug in the power adapter.
The battery low icon or the bAt Lo message are displayed on the screen	The battery is low.	The battery is low and needs to be recharged.
out	The blood pressure measurement is out of the measurement range.	Re-measure and follow the instructions.
E1	The cuff is not tied tightly or inflated abnormally.	Re-adjust the armband, make it moderately tight, and measure again.
E2	Measurement errors are caused by hand movement, conversation or weak arterial pulse during measurement.	Relax and re-measure.
E3	No pulse signal was detected during the measurement.	Loosen the clothes on the arm and measure again.
E4	Blood pressure measurement failed.	Re-measure after relaxing and adjusting.
EE XX	Calibration error (XX can be some digital symbols such as 01, 02, etc., and similar situations are all calibration errors)	Re-measure. If the problem remains unresolved, please contact the dealer or our customer service department for further assistance. Please refer to the instructions on the warranty card for the method of return or repair.

Maintenance

Cleaning/disinfection process:

Step 1: Make sure to switch off and unplug the device prior to cleaning.
Step 2: Use a soft cloth wetted with soapy water to clean the cuff first, and then use a soft cloth wetted with clear water to remove residual soap until there is no visible residual contaminants. To disinfect the cuff use a soft cloth wetted with 70% isopropanol for about 3 minutes in case of domestic use and for about 10 minutes in case of professional use. Attention shall be paid to avoid liquid penetration into the cuff.

Step 3: Use a dry soft cloth to wipe the cuff, in order to remove residual moisture.

Step 4: Dry the cuff at a well-ventilated place after cleaning.

Suggestion:

Frequency of cleaning and disinfection: for single patient multiple use, it's recommended to clean the device surface once a month or whenever it's necessary; for multiple patient multiple use, it's recommended to clean the device every time before and after usage. Maintenance procedures shall be taken as per instruction.

Calibration and service

The accuracy of this monitor has been carefully tested and is has been designed for a long life. Indeed, the monitor can maintain the safety and performance characteristics for a minimum of 10,000 measurements or five years of normal use.
It is generally recommended to check the unit every two years to ensure correct functioning and accuracy. Please consult the contact information in this user manual to get in touch with your local authorized Service Center.

WARNING

- Report any serious accident occurring in connection with the device to your competent authority and to the manufacturer.
- Do not confuse self-monitoring with self-diagnosis. By using this device you can monitor your blood pressure, but it does not replace your doctor. Begin or end medical treatment based solely on medical advice. Never change a prescribed medication without consulting your physician.
- If you are taking medication, consult your physician to determine the proper time to measure your blood pressure.

- Store your device, cuff and adapter in a clean and dry place, protect it against extreme moisture, heat, lint, dust and direct sunlight. Never place any heavy objects on it.
- Keep out of reach of children, people with reduced cognitive abilities or pets: the device contains small parts which can cause suffocation or internal lesions if swallowed. The USB cable can create strangulation risk.
- If the arm cuff pressure reaches 300 mmHg it deflates automatically. If the arm cuff does not deflate when pressure reaches 300 mmHg or if you feel discomfort during the measurement take it off your arm and press the START/STOP button to stop inflation.
- If Irregular Heartbeat (IHB) brought by common arrhythmias is detected in the procedure of blood pressure measurement, a signal of will be displayed. Under this condition, the device can keep function, but the results may not be accurate, it's suggested that you consult with your physician for accurate assessment.
- This unit is not suitable for continuous monitoring during medical emergencies or operations. Otherwise, the patient's arm and fingers will become anesthetic, swollen and even purple due to a lack of blood.
- This device cannot be used with HF surgical equipment at the same time.
- This device is not intended to be used during patient transportation outside a healthcare facility.
- It is not intended for use on extremities other than the arm or for functions other than obtaining a blood pressure measurement.
- To avoid measurement errors, please avoid the condition of strong electromagnetic field radiated interference signal or electrical fast transient/burst signal (see EMC tables at the end of these IFU for further details).
- The user must check that the equipment functions safely and see that it is in proper working condition before being used.
- Manufacturer will make available on request circuit diagrams, component parts list etc. only to qualified person.
- Too frequent measurements can cause injury to the patient due to blood flow interference. Please check (by observation of the limb concerned) that operation of OneRAPID does not result in prolonged impairment of the circulation of the blood of the patient.
- Please use the device under the environment which was provided in the instruction manual. If this is not complied with, correct operation and life of the device may be compromised.
- During using, the patient will be in contact with the cuff. The materials of the cuff have been tested and found to comply with requirements of EN ISO 10993-5:2009, EN ISO 10993-10:2021 and EN ISO 10993-23:2021. It will not cause any potential sensitization or irritation reaction.
- The application of the cuff over a wound can cause further injury. Application is not recommended.
- The application of the cuff and its pressurization on any limb where intravascular access or therapy, or an arterio-venous (A-V) shunt, is present because of temporary interference to blood flow and could result in injury to the patient.
- The application of the cuff and its pressurization on the arm on the side of a mastectomy can cause injury or give not accurate measurement.
- Do not inflate the cuff on the same limb which other monitoring ME equipment is applied around simultaneously, because this could cause temporary loss of function of the simultaneously-used monitoring ME equipment.
- The device is a non-AP/APG device and not suitable for use in the presence of anaesthetic mixtures that are inflammable with air, oxygen or nitrous oxide.
- The operator shall not touch DC output jack of AC/DC adapter and the patient simultaneously.
- Please use accessories and detachable parts specified/authorized by the manufacturer. Otherwise, it may cause damage to the unit or danger to the user/patients.
- We recommend checking its performance every two years or after it has been repaired. The device must be checked again in case of damage due to impact (such as falls) or exposure to liquids and/or extreme temperatures (hot/cold) or extreme humidity variations
- Self-use, without supervision, in public place is not allowed.
- Please dispose of accessories, detachable parts, and the ME EQUIPMENT according to the local guidelines.
- In case of problems concerning performances, settings, maintenance or use of the device, contact your local distributor or the Italian toll-free number 800 900 080. Warning! Do not service and/or repair the appliance while it is being used! Any maintenance/service activity/operation must be done solely by a Pic Service Centre. No modification of the equipment is allowed.
- When the device is stored at the minimum/maximum storage temperature it takes 4h to warm/cool until the device is ready for its intended use.
- The maximum temperature that the applied part can be achieved is 41.9 °C while the environmental temperature is 40 °C and the contact time with the cuff by patient should be less than 10 minutes.

	Follow the Instruction for Use		BF-type Applied Part
	Symbol of compliance with Regulation (EU) 2017/745. Compliance to RED directive 2014/53/EU		Recycling (cardboard)
	Manufacturer		Recycling (plastic)
	Direct Current		Production date
	Latex free		Recycling (plastic)
	Atmospheric pressure limitation		Recycling (plastic)
	Humidity limitation		Recyclable packaging
	Temperature limitations		Triman
	Warning		Der Grüne Punkt

	Batch number		Medical Device
	The Bluetooth® combination mark, Bluetooth® version v. 5.0		Catalogue number
	Unique Device Identifier		Serial number
	3,6V -1000mAh Li-ion battery		AC/DC Adapter
	The product must be delivered to a separate waste collection centre for electrical and electronic devices or returned to the retailer when purchasing a new equivalent device.		
	This product meets the basic safety and essential performance requirements indicated in the IP22 conditioning test (protection against solid foreign objects of 12.5 mm in diameter and greater vertically falling water drops when the enclosure is tilted up to 15°).		

SPECIFICATIONS

Power supply:	3,6V -1000mAh Li-ion battery. External USB connector for battery charging only (with all device's functions disabled).
Display:	V.A. 72mm x 22mm
Measurement mode:	Oscillometric inflation -type measurement mode
Measurement range:	Rated cuff pressure: 0mmHg~299mmHg
Measurement pressure:	SYS: 60-230mmHg, DIA: 40-130mmHg
Pulse value:	40-199 beat/minute
Accuracy:	Pressure: (5°C-40°C) ± 3mmHg Pulse value: ±5%
Normal working condition:	Temperatures: +5°C+40°C Relative humidity: 15%~90%, Atmospheric pressure: 700hPa to 1060hPa
Storage & transportation condition:	Temperatures: -20°C~+60°C Relative humidity: ≤93%, Atmospheric pressure: 500 hPa to 1060 hPa
Arm circumference:	22cm~42cm
Net Weight:	Approx 204 g
External dimensions:	123x44x22mm
Mode of operation:	Continuous operation
Useful life:	10,000 Measurements
Degree of protection:	Type BF applied part (cuff)
Protection against penetration of solids and liquids:	IP22 - Protected against penetration by solid foreign bodies of 12.5 mm or more and protection against penetration of liquids (dripping water when tilted 15°) A.01.01 AC/DC Adapter: 100-240V~ 50/60Hz / 5V 1A AC/DC INPUT 100-240V ~ 50/60 Hz 0,2 A max OUTPUT 5V 1000 mA. Adaptor REF 02010309100000 (EU) or REF 02010308100000 (UK) please refer to customer service
Software version:	
Accessories:	

THIS PRODUCT COMPLIES WITH THE DIRECTIVE 2012/19/EU.

The crossed bin symbol on the appliance indicates that the product, at the end of its life, must be disposed of separately from domestic waste, either by taking it to a separate waste disposal site for electric and electronic appliances or by returning it to your dealer when you buy another similar appliance. The user is responsible for taking the appliance to a special waste disposal site at the end of its life. If the disused appliance is collected correctly as separate waste, it can be recycled, treated and disposed of ecologically; this avoids a negative impact on both the environment and health, and contributes towards the recycling of the product's materials. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

THIS PRODUCT COMPLIES WITH EU DIRECTIVE 2006/66/EC

The crossed bin symbol on the batteries indicates that, at the end of their life, they must be disposed of separately from domestic waste, either by taking them to a separate waste disposal site for batteries or by returning them to your dealer when you buy similar rechargeable or non-rechargeable batteries. The chemical symbols Hg, Cd, Pb, printed under the crossed bin symbol, indicate the type of substance contained in the batteries: Hg=Mercury, Cd=Cadmium, Pb=Lead. The user is responsible for bringing batteries, at the end of their life, to the appropriate collection facilities in order to facilitate treatment and recycling. The correct disposal of your old toy will help prevent potential negative consequences for the environment and human health and favour the recycling of the materials of which the product is made. Illegal disposal of the product by the user will damage the environment and human health. For further information regarding the waste disposal services available, contact your local waste disposal agency or the shop where you bought the appliance.

DECLARATION OF EC COMPLIANCE:

PIKDARE S.p.A. hereby declares that this appliance, OneRAPID complies with all the essential requirements and other provisions set forth by the EU Directives: 2014/53/EU (RED), 2011/65/EU (RoHS) and EU 745/2017 (MDR). The full text of the EU declaration of conformity is available at the following internet address: <http://web.picsolution.com/pichealthstation/en/onerapid/en.html>. In accordance with European Commission Decision No. 2000/299/EC dated 06/04/2000, the frequency band used by this product has been harmonised for all EU countries and this is therefore a Class 1 product which can be used in all European Union countries.

Pouch cleaning symbols

- Wash by hand in cold water,
- Do not bleach,

- Do not tumble dry,
- Line drying,
- Do not iron,
- Do not dry clean.

PIKDARE S.p.A.
Via Saldarini Catelli, 10
22070 - Casnate con Bernate (CO) - Italy
www.picsolution.com

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EMC Guidance

- This product needs special precautions regarding EMC and needs to be installed and put into service according to the EMC information provided, and this unit can be affected by portable and mobile RF communications equipment.
- Do not use a mobile phone or other devices that emit electromagnetic fields, near the unit. This may result in incorrect operation of the unit.
- Caution: This unit has been thoroughly tested and inspected to assure proper performance and operation!
- Caution: This machine should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, this machine should be observed to verify normal operation in the configuration in which it will be used.

Guidance and manufacturer's declaration - electromagnetic emissions		
The OneRAPID (REF) 02010394000000 is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class [B]	
Harmonic emissions IEC 61000-3-2	Class A	The OneRAPID (REF) 02010394000000 is suitable for use in all establishments, other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration - electromagnetic immunity			
The OneRAPID (REF) 02010394000000 is intended for use in the electromagnetic environment specified below. The customer or the user of the OneRAPID (REF) 02010394000000 should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±15 kV air	±8 kV contact ±15 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	Power supply lines: ±2 kV Input/output lines: ±1 kV	Power supply lines: ±2 kV	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	line(s) to line(s): ±1 kV line(s) to earth: ±2 kV 100 kHz repetition frequency	line(s) to line(s): ±1 kV 100 kHz repetition frequency	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	0% UT; 0,5 cycle A 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase; at 0° 0% UT; 300 cycle	0% UT; 0,5 cycle A 0°, 45°, 90°, 135°, 180°, 225°, 270° and 315° 0% UT; 1 cycle and 70% UT; 25/30 cycles Single phase; at 0° 0% UT; 300 cycle	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50/60Hz) magnetic field IEC 61000-4-8	30 A/m 50Hz/60Hz	30 A/m 50Hz/60Hz	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Guidance and manufacturer's declaration - electromagnetic immunity			
The OneRAPID (REF) 02010394000000 is intended for use in the electromagnetic environment specified below. The customer or the user of the OneRAPID (REF) 02010394000000 should assure that it is used in such an environment.			
Immunity Test	IEC 60601 Test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	150 kHz to 80 MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	150 kHz to 80 MHz: 3Vrms 6Vrms (in ISM and amateur radio bands) 80% Am at 1kHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables, than the recommended separation distance calculated from the equation appropriate for the frequency of the transmitter. Recommended separation distances: $d = 0,35 \sqrt{P}$ 80 MHz to 800 MHz: $d = 1,2 \sqrt{P}$ 800 MHz to 2,7 GHz: $d = 2,3 \sqrt{P}$ Where p is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. ¹ Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	10V/m, 80% Am at 1kHz	10V/m, 80% Am at 1kHz	

- NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

¹Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the device is used exceeds the applicable RF compliance level above, the device should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the device.
²Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

Recommended separation distances between portable and mobile RF communications equipment and the device.			
The OneRAPID (REF) 02010394000000 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the OneRAPID (REF) 02010394000000 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the device as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=3,5 \sqrt{P}$	80MHz to 800MHz $d=1,2 \sqrt{P}$	800MHz to 2,7GHz $d=2,3 \sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.
NOTE 1: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.
NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration - electromagnetic immunity						
The OneRAPID (REF) 02010394000000 is intended for use in the electromagnetic environment specified below. The customer or the user of the OneRAPID (REF) 02010394000000, should assure that it is used in such an environment.						
Test frequency (MHz)	Band *) (MHz)	Service *)	Modulation *)	Modulation *) (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation *) 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	FM *) ±5 kHz deviation 1 kHz sine	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation *) 217 Hz	0,2	0,3	9
745						
780						
810	800 - 960	GSM 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation *) 18 Hz	2	0,3	28
870						
930						
1 720	1 700 - 1 990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25; UMTS	Pulse modulation *) 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2 400 - 2 570	Bluetooth*, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation *) 217 Hz	2	0,3	28
5 240						
5 500						
5 785	5 100 - 5 800	WLAN 802.11 a/n	Pulse modulation *) 217 Hz	0,2	0,3	9

NOTE: If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.

- For some services, only the uplink frequencies are included.
- The carrier shall be modulated using a 50% duty cycle square wave signal.
- As an alternative to FM modulation, 50% pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.

The MANUFACTURER should consider reducing the minimum separation distance, based on RISK MANAGEMENT, and using higher IMMUNITY TEST LEVELS that are appropriate for the reduced minimum separation distance. Minimum separation distances for higher IMMUNITY TEST LEVELS shall be calculated using the following equation:

$$E = \frac{P}{d} \sqrt{P}$$

Where P is the maximum power in W, d is the minimum separation distance in m, and E is the IMMUNITY TEST LEVEL in V/m.

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